Zenthera Zapper v1.2 User Guide

Firmware version 1.2#20260701

Dated : 01 Jul 2026

Table of Contents

Discover the Power of Zenthera Zapper 5
Safety Guidelines 6
Battery Care Instructions 7
- Airline Travel Precautions 7
- Battery Status Indicators 7
- Charging Instructions 8
- Pro Tips 8

User Operating Manual 8
- Welcome to Zenthera Zapper 8

Quick Start Guide 9
- 1. Power On 9
- 2. Select Mode 10
- 3. Start Session 10
- 4. During Session 10

Display Explained 10
Using Carbon Electrodes (Pads / Straps) 11
Core Modes Summary 12
- Daily Essentials (Use Electrodes) 12
- Specialized Health (Use Electrodes) 12
- Additional Health Modes (Use Electrodes) 13

Sound Healing Modes (Use External Speaker with AUX Cable - No Electrodes) 14
Special Modes 14
03 Sweep Mode (60 min) 14
26 ParasiteSweep (30 min) 14
31 Manual Mode Hz (7 min) 14
51 Manual Mode KHz (7 min) 15
Button Functions by State 15
Complete Session Flow Example 15
Daily Usage Recommendations 16
Beginner (30 min/day) 16
Regular / Maintenance (45–120 min/day) 16
Intensive Protocol 16
External Speaker with AUX Cable and Volume Control Usage (Sound Healing Modes 19-22, 27, 28) 17
Selecting the external Speakers for Sound Healing Modes 17
Must-have specifications 17
Recommended (best results) specs 18
What to avoid 18
Troubleshooting 18
Frequency Table by Mode 19
Detailed Sweep & Manual Mode Explanations 25
Mode 03: Sweep Mode 25
Mode 26: Parasite Sweep 25
Mode 31: Manual Mode 26
Special Modes (Non-Medical / Sound Healing) 26
01 Program Mode 28
Detailed Illness Associations and Frequency Explanations 28
Step 1: 2128 Hz (Dr. Piyush Saxena's Primary Frequency) 28
Step 2: 10 kHz (Systemic Antimicrobial) 29

Step 3: 30 kHz (Dr. Hulda Clark Universal Killer) 29
Step 4: 48 kHz (Deep Systemic Reset) 29
Usage Recommendations for Program Mode 30
30 Factory Reset 30
Women's Health 30
33 UTI Kidney 31
34 Skin Wound 31
35 Anemia Care 31
36 Cervix HPV 31
37 Menstrual Care 31
38 Vaginitis 31
39 PID Adnexitis 31
40 Urinary Incont 31
41 STIs RTIs 32
42 Ovarian Cysts 32
43 Menopause Help 32
Other Specialized Health Modes 32
44 Thyroid Care 32
45 Bone Health 32
46 GallBldr Pain 32
47 Migraine Relief 32
48 Sinus/Vertigo Focus 32
49 Eyes Support 33
50 Liver Support 33
Product Warranty (Hardware) 48
Limited Hardware Warranty 48
Warranty Claim Process: 49
Disclaimer of Medical Warranty 49

Limitation of Liability 50
Safety Disclaimer (Critical) 50
Software/Firmware Disclaimer 51
No Endorsement Disclaimer 52

Whats New in Firmware version 1.2#20260701

- Splited “31 Manual Mode” in to “31 Manual Mode Hz” and “51 Manual Mode KHz” for easy custom frequency setting.
- Added few more frequencies in “14 Parasite Care” and “50 Liver Support”.
- Release date 01 July 2026
- Update Type: Factory Update

Discover the Power of Zenthera Zapper



Zenthera Zapper brings ancient healing wisdom into modern technology with **50 precision-programmed frequency modes** that target everything from common infections and chronic pain to energy balancing and pet wellness. Inspired by Dr. Royal Rife's frequency therapy, Dr. Hulda Clark's protocols, and Dr. Piyush Saxena's cleansing sequences, this compact device delivers targeted electrical frequencies (Rife modes) or soothing sound vibrations (Sound Healing modes) through simple electrode or External Speaker with AUX Cable connections. Whether you're seeking daily detoxification, viral support, joint pain relief, or deep relaxation, Zenthera Zapper automates complete treatment protocols with crystal-clear LCD feedback showing real-time frequency, countdown timer, and progress.

Zenthera Zapper is safe for Humens / Pets and Food, it runs on 5v to 6v DC current. Which will not feel until you use it on thin skin like stomach.

No guesswork, no complex setup – just press MODE and heal. With intuitive 3-button operation, auto-advancing 7-minute sessions, and smart features like frequency skipping (UP/DOWN) and session memory, Zenthera Zapper fits seamlessly into your wellness routine. Run the **28-minute** Program Mode daily for systemic cleansing, choose from a wide range of electrode-based health protocols (roughly 21–240 minutes, depending on the mode), or relax with Sound Healing modes (about 5 -112 minutes) featuring Chakra tones, Brainwave entrainment, and Solfeggio-style frequencies. Safe, non-invasive, and complementary to medical care, Zenthera Zapper empowers you to take charge of your natural healing journey – one frequency at a time.

Safety Guidelines

- **Use as directed** - complementary to medical care
- **3 / 7-minute sessions** prevent fatigue
- **Rest 20+ minutes** between modes
- **Hydrate well** during use
- **NOT for**: Pacemakers / pregnancy (consult doctor)
  Stop if: Dizziness, discomfort occurs
- **Medical Disclaimer**: Zenthera Zapper provides experimental frequencies. This device provides experimental wellness frequencies; **no guaranteed results**; and not medical device.
- **Implants/medical devices:** Do not use if you have a pacemaker/ICD, implanted neurostimulator, insulin pump, cochlear implant, or any implanted electronic/metal device unless your doctor approves.
- **Heart/neck placement:** Never place electrodes across the chest/heart (left hand to right hand) or on the front of the neck.
- **Children & supervision:** Not for children without adult supervision/doctor guidance.
- **Skin safety:** Do not use on broken skin/open wounds with electrodes; stop if burning/tingling/pain occurs; clean electrodes after use.
- **Driving/sleep warning for sound modes:** Do not use Sleep/Brainwave sound healing modes while driving / operating machinery.
- **No replacement for treatment:** Strengthen to: Do not delay emergency care (chest pain, breathlessness, high fever, heavy bleeding, severe abdominal pain).
- **If symptoms worsen:** Stop using and consult a doctor if symptoms worsen or persist >48 hours.
- **Charging while using:** Do not use the device while charging.

- **Heat/fire:** Do not leave charging unattended, do not charge near pillow/blanket, stop charging if battery/device gets hot or swollen.

Battery Care Instructions

Power Source: 1000 mAh internal Li-ion battery

Airline Travel Precautions

1. **Check airline policy** before flying. Most allow removable batteries <100Wh (Watt-hours).
2. **Remove battery for security**:
    - Unscrew 4 screws (keep loose, don't tighten)
    - Disconnect battery clip, remove battery
    - Place **Zapper body + loose screws** in one transparent zip bag
    - Place **Battery** in separate transparent zip bag
3. **No screwdriver** in carry-on/handbag
4. **Security demo**: Show how battery installs if asked

Battery Status Indicators

- LCD contrast low → Battery ~10%
- LCD barely visible → Battery ~5%
- LCD completely off → Battery 0% (dead)

Note: The zapper shows "**Batt Left: XX%**" at boot and has a critical low-battery lockout. when battery is repeatedly detected below the required power, the device beeps and displays "**Battery LOW / CHARGE NOW**". This is the time to charge the battery.

Charging Instructions

1. **Power OFF** device completely
2. Connect **9V 0.5A-1A AC-DC adapter** to power jack (use a good-quality 9V adapter (stable output). Avoid cheap unbranded adapters.)
3. Charge for **~2.0 hours maximum** - **avoid overcharging**
4. Disconnect charger when done


Pro Tips

- Charge when LCD contrast fades
- Keep spare battery for long sessions
- Store at 50-60% charge for longevity
- Charge every 2-3 months if unused

User Operating Manual

Welcome to Zenthera Zapper

Zenthera Zapper delivers precise healing frequencies through electrodes (health modes) or sound vibrations (Sound Healing modes). With **50 specialized modes**, it covers everything from infections and pain relief to chakra balancing and pet care. Each mode runs automatically with a clear LCD display showing frequency, time remaining, and progress.



Simple 3-Button Operation:

- **MODE button**: Select modes / Start-Stop sessions
- **UP button**: Previous mode / Next frequency
- **DOWN button**: Next mode / Previous frequency

Quick Start Guide

1. Power On

"`Zenthera Zapper`" and its firmware details ← Shows for 2 seconds

"`DO NOT USE WITH PACEMKR/PREGNNCY`" ← Warning will be displayed for 2.5 seconds.

"`Batt Left: ~XX%`" ← Battery screen for 1.5 seconds

"`Select Mode: 01 Program Mode`" ← Automatically shows.

Press **UP** or **DOWN** to change modes, **MODE** to start session

2. Select Mode

1. Press **MODE** → Shows menu: Select Mode: "`01 Program Mode`"
2. Use **UP/DOWN** to browse different modes
3. Last used mode remembers position

3. Start Session

Press **MODE** again → Session begins automatically:

`Vir: 000.728 kHz` ← Mode + Current Frequency

`Tim:07:23 Rem:13` ← Time remaining + Frequencies left

4. During Session

- **UP**: Skip to next frequency (resets 7-min timer)
- **DOWN**: Go back to previous frequency
- **MODE**: Stop session (long beep)

Display Explained

`Com: 000.727 kHz` ← "Com" = Common Toxins mode, 727 Hz

`Tim:06:45 Rem:12` ← 6 min 45 sec left, 12 frequencies remaining

`Chk: 432 Hz OM` ← Chakra mode showing "OM" label

`Tim:07:12 Rem:0` ← Last frequency in set

- **Tim**: Remaining time for current frequency
- **Rem**: Frequencies left in complete set (most modes)



Using Carbon Electrodes (Pads / Straps)

- Use the supplied **carbon rubber electrodes** for all Health Modes that require electrodes (01-18, 02, 23-26, 31-47). The Sound Healing Modes 19–22, 27, 28 use an external speaker with AUX cable, not electrodes.
- Each pad has a **black carbon side** and a **blue plastic (insulated) side**. The black carbon side must touch the skin; the blue plastic side always faces outward.
- Apply a thin, even layer of conductive gel only on the **black carbon side** before placing it on the skin. Do not apply gel on the blue plastic side.
- Skin should be clean and dry (no oil, lotion, or broken skin). Do not place pads on cuts, rashes, infected or irritated areas. If you feel strong burning, sharp pain, or unusual discomfort under any pad, stop the session immediately and adjust pad position, gel, or end the session.
- For general whole-body support (Program, Common Toxins, Parasite Care, Detox Body, Immune Boost, etc.), place one pad on each wrist or each ankle so that the current passes through the body between the two pads.
- For local support (knee, shoulder, back, specific joint or area), place both pads so that the line between them “passes through” the problem area (for example, either side of the knee, above and below a painful spot on the spine, or front and back of a shoulder).
- Do not place electrodes **across the chest** (for example, one pad on the left chest and one on the right chest), on the front of the neck, or on the head/face.

- Use only on intact skin. Do not use if you have a pacemaker or are pregnant; consult your doctor before use for any serious medical Condition. Stop use if you experience dizziness, chest discomfort, breathing difficulty, or any unusual symptoms.
- After each session, gently remove the pads, wipe off excess gel from the skin and pad surface with a slightly damp cloth, squeeze out extra water from the pads if you rinsed them, and let the pads air-dry before storing. Ensure the device output jacks and cables stay dry at all times.

Core Modes Summary

Daily Essentials (Use Electrodes)

Mode	Purpose	Total Time	Best For
01 Program	Dr. Piyush Saxena Protocol	**28 min**	Daily cleanse
04 Common	Cold, fever, infection	**111 min**	General health
05 Viral	Flu, COVID, dengue	**98 min**	Viral support
23 Detox	Body cleansing	**35 min**	Weekly detox

Specialized Health (Use Electrodes)

Mode	Target	Frequencies	Total Time
06 ENT	Ear, sinus, throat	12	84 min
07 Dental	Toothache, gums	13	91 min
08 Cardiac	Heart, circulation	11	77 min
13 Arthritis	Joint pain, back	12	84 min
14 Parasite	Worms, malaria	23	161 min
29 Rife Classic	Classic Frequencies of Rife Protocol	11	77 min
02 Core Support	Complete Core Frequency Set	121	242 min

Additional Health Modes (Use Electrodes)

Mode	Target	Frequencies	Total Time
32 Pemphigus Vul	Autoimmune skin blistering support	14	98 min
33 UTI Kidney	UTI / E. coli urinary support	36	108 min
34 Skin Wound	Skin/wound infection support	10	70 min
35 Anemia Care	Anaemia / blood support set	43	129 min
36 Cervix HPV	Cervical HPV support	16	112 min
37 Menstrual Care	Menstrual/endometriosis/fibroid support	18	126 min
38 Vaginitis	Vaginitis/vaginosis/Candida support	22	154 min
39 PID Adnexitis	PID/adnexitis support	7	49 min
40 Urinary Incont	Urinary control/bladder support	11	77 min
41 STIs RTIs	STI/RTI support set	15	105 min
42 Ovarian Cysts	Ovarian cyst support	3	21 min
43 Menopause Help	Menopause support	16	112 min
44 Thyroid Care	Thyroid support	11	77 min
45 Bone Health	Bone regeneration support	16	112 min
46 GallBldr Pain	Gallbladder support	41	123 min
47 Migraine Relief	Migraine / Headache relief	5	35 min
48 Sinus/Vertigo	Sinus, Sinusitis, Rhinitis, Tinnitus, Vertigo Support	93	186 min
49 Eyes Support	Cataracts, Glaucoma, Macular Degeneration, Eye Strain, Dry Eyes, Diabetic Retinopathy, Conjunctivitis, Blurred vision	63	186 min
50 Liver Support	Detox Of Anesthesia, Fatty Liver, Hangover, Jaundice, Lipoma, Hepatitis (A/B/C/Chronic/General), Listerioses, Liver Cirrhosis, Liver	120	240 min

Mode	Target	Frequencies	Total Time
	Diseases, Liver Enlargement, Parasites Flukes Liver, Pesticide Detox, Vitiligo, Yellow Fever		

Sound Healing Modes (Use External Speaker with AUX Cable - No Electrodes)

Mode	Purpose	Total Time
19 7-Chakras	Energy balancing	42 min
20 Beej Mantra	Subtle energy	49 min
21 Brain Wave	Sleep, focus	98 min
22 Healing	Pain relief, DNA repair	112 min
27 Pets & Food	Pets/water energizing	5-30 min
28 Sleep & Relax	Helps to Sleep Stress relief	55 min

Note: Many frequencies (2-45 Hz) are too low to hear but work through vibration.

Special Modes

03 Sweep Mode (60 min)

- **Automatic scan**: 2 Hz → 500 kHz
- Covers full healing spectrum

26 ParasiteSweep (30 min)

- **Specialized scan**: 20 Hz → 50 kHz
- Targets parasite frequency range

31 Manual Mode Hz (7 min)

- This mode is designed for the range of 2 Hz to 999 Hz with the default 50 Hz allows you to set any frequency in between for 7 min.

- Start: 50 Hz
- **UP**: Increases by 25 Hz at the max 999 Hz, after 999 Hz back to 50 Hz
- **DOWN**: Decreases by 1 Hz
- Perfect for custom frequencies, if you want to reset the frequency back to 50 Hz, then follow the **Factory Rest** process.

51 Manual Mode KHz (7 min)

- This mode is designed for the range of 1 KHz to 500 KHz with the default 30 KHz allows you to set any frequency in between for 7 min.
- Start: 30 KHz
- **UP**: Increases by 25 KHz at the max 500 KHz, after 500 KHz back to 30 KHz
- **DOWN**: Decreases by 1 KHz
- Perfect for custom frequencies, if you want to reset the frequency back to 30 KHz, then follow the **Factory Rest** process.

Button Functions by State

State	MODE	UP	DOWN
Idle	Enter Menu	-	-
Menu	Start Session	Previous Mode	Next Mode
Running	STOP (long beep)	Next Freq	Previous Freq

Complete Session Flow Example

1. Power ON → "Press MODE to Start"
2. MODE → "Select Mode: 01 Program Mode"
3. UP/DOWN → Find desired mode

4. MODE → "Prg: 002.128 kHz" / "Tim:07:00 Rem:3"
5. Waits 7 min → Auto advances: "Prg: 10.000 kHz"
6. UP anytime → Skip to next frequency immediately
7. After 4th frequency → "Program Done!" → Back to Start

Daily Usage Recommendations

Beginner (30 min/day)

Day 1: 01 Program Mode (28 min)

Day 2: 04 Common Toxins (run only first 4 steps ≈ 28 min; use **UP** to skip to the next frequency)

Day 3: 23 Detox Body (35 min)

Note: If you feel discomfort, press MODE to stop and resume next day

Regular / Maintenance (45–120 min/day)

Option A (short): 01 Program (28 min) + one Health Mode for 15–60 min (use **UP** to shorten long sessions).

Option B (full): 01 Program (28 min) + one full Health Mode session (typical 60–120 min depending on mode).

Note: Long modes are optional to finish in one sitting; you can stop anytime with MODE and continue later.

Intensive Protocol

Morning: 01 Program (28 min) + 14 Parasite Care (shortened using **UP**, e.g., run 4–8 steps first; do not try to complete the whole mode in one go initially).

Evening: 04 Common Toxins (shortened if needed using **UP**) + any preferred Sound Healing mode (19–22, 28) using external speaker (no electrodes).

Note: Rest 20 minutes between long sessions and hydrate well.

External Speaker with AUX Cable and Volume Control Usage (Sound Healing Modes 19-22, 27, 28)

Disconnect electrodes completely. Connect External Speaker with AUX Cable/speaker. These gentle frequencies (2-963 Hz) work through calming sound vibrations. Many are subsonic (too low to hear) but powerful for relaxation, chakra balance, and pet calming. Use headphones for personal sessions or speakers for room therapy.

Mode	Volume	Speaker Type	Best Setup
19-20 Chakras/Beej	**Low-Medium**	Headphones	Personal energy work
21 Brain Wave	Medium	Stereo speakers	Binaural entrainment
22 Healing	Medium-Loud	Room speakers	Solfeggio immersion
27 Pets & Food	Loud	Bowl next to speaker	Vibrational transfer
28 Sleep & Relax	Low-Medium	Headphones or Stereo speakers near the bed	Restful slumber, stress relief

Selecting the external Speakers for Sound Healing Modes

External speakers for Zenthera Zapper's Sound Healing modes (19–22, 27& 28) should be chosen to reproduce 2 Hz – 1 kHz reliably at comfortable volume, using the External Speaker with AUX Cable (no electrodes).



Must-have specifications

- **Input**: 3.5 mm AUX/line-in (stereo preferred).
- **Speaker type**: Active (powered) speakers or any speaker connected to a small amplifier (phone-passive speakers alone won't work).
- **Frequency response** (realistic): At least 40 Hz - 10 kHz on paper; lower is better, but most consumer speakers won't truly reproduce 2–20 Hz (those will be felt as vibration more than heard).

- **Low noise floor**: minimal hiss at low volume, since these modes are used at low–medium volume and often during relaxation/sleep.
- **Volume Control**: The Speaker must have volume control.

Recommended (best results) specs

- **Stereo playback**: For Brain Wave (Mode 21) and Sleep/Relax-style entrainment, use stereo speakers (L/R) or stereo headphones, because the intended setup for brainwave entrainment is stereo.
- **Bass capability**: a speaker system that can reach ~40 Hz (or a 2.0 + subwoofer system reaching ~30–35 Hz) helps you feel the low-frequency components more clearly.
- **Power**: roughly 10–50 W RMS total is enough for a bedroom/room at low to medium volume without distortion.

What to avoid

- Bluetooth-only speakers with no AUX input, because the document's recommended connection method is the AUX cable.
- Very small/mono speakers if you want the full "Brain Wave" style effect; stereo is preferred for those modes.

Troubleshooting

- No display? → Battery Drained, turn the on/off button to OFF and connect the 9v AC-DC Adapter. Wait for 90 min.
- No sound/beeps? → Device working normally
- Session won't start? → Press MODE twice
- Frequency stuck? → Use UP or DOWN button

Support: Last mode auto-saves. Power cycle returns to your favourite mode.

Frequency Table by Mode

Note: Most modes run 7 minutes per step; some modes run **3 minutes** per step (e.g., Common Toxins, Parasite Care, Cancer Support, Rife Classics, UTI+Kidney, Anemia, Gallbladder, Eyes Support), and some run **2 minutes** per step (Core Support, Sinus/Vertigo).

Mode	Illnesses Targeted	Rationale / Explanation	Total Time in minutes	Auto - Advance	Short Form
01 Program Mode	General Protocol (Dr. Piyush Saxena)	High-frequency protocol steps; 2128 common in systemic cleansing, 30kHz from Dr. Clark/Hulda Clark protocols	28	Yes	Prg
02 Core Support	Broad-spectrum helper set	Complete backbone frequency set. 2 Minutes per frequency, total 121 frequencies.	242	Yes	CoS
03 Sweep Mode	Broad-Spectrum Pathogen Scan	Full spectrum sweep technique used in Rife experimentation to hit unknown pathogens across entire frequency range	60	Continuous	Swp
04 Common Toxins	Cough, Cold, Fever, Headache, Sore Throat, Allergy, Asthma, Bronchitis, Sinusitis, Diarrhea, Nausea, Fungus, Boils/Pimples (Skin Infection), General Infection	37 frequencies running, 3 minutes each frequency.	111	Yes	Com
05 Viral Relief	COVID, Dengue, Flu, Malaria, Hepatitis, Typhoid, Gastroenteritis	Classic antiviral trio + high harmonics for lung/viral load	98	Yes	Vir
06 ENT Care	Ear Infection, Tinnitus, Sinus, Sore Throat, Tonsillitis, Vertigo	ENT infection core and tinnitus/hearing support	196	Yes	ENT

Mode	Illnesses Targeted	Rationale / Explanation	Total Time in minutes	Auto - Advance	Short Form
07 Dental Care	Toothache, Gums, Abscess, Oral Thrush, Jaw Pain	Dental abscess protocol (465-880 core); 10000 for chronic infection	126	Yes	Den
08 Cardiac Care	Arrhythmia, Angina, Blood Pressure, Circulation	Arrhythmia-specific (225/510); vascular support harmonics	77	Yes	Car
09 Metabolic Care	Diabetes, Obesity, Thyroid, Kidney, Liver	Infection + liver/kidney detox; terrain support for metabolic stress	98	Yes	Met
10 Organ Health (General)	Lungs, Stomach, Brain, Eyes, Skin, Pancreas, COPD	Multi-organ infection + detox; lung/brain inflammation. 3 min each frequency	165	Yes	Org
11 Reproductive Organs	PCOS, Fibroids, UTI, HPV, Prostate, Infertility	Pelvic infection core; structural support for fibroids/cysts	112	Yes	Rep
12 Autoimmune	Lupus, Psoriasis, Crohn's, Vitiligo, Celiac	Immune modulation + infection terrain; not disease-specific	84	Yes	AIm
13 Arthritis Care	Joint/Back Pain, Sciatica, Fibromyalgia, Gout, OA/RA	Pain relief + arthritis core	84	Yes	Art
14 Parasite Care	Worms, Giardia, Amoeba, Tapeworm, Malaria Parasite	Low fundamentals + high parasite harmonics (tapeworm range) – 3 min each frequency	195	Yes	Par
15 Cancer Care	Breast, Lung, Prostate, Colon, Lymphoma Support	Cancer terrain + infection/pain support - 3 min each frequency	84	Yes	Can
16 Ulcer Care	Peptic Ulcers, Stomach Ulcers, Fistula support (Dental fistula, Tracheoesophageal fistula, Fistula general), General Ulcer	H. pylori/infection core + stomach lining support; includes fistula/dental fistula support frequencies	81	Yes	Ulc
17 Piles & Gut	Haemorrhoids, Gut Inflammation, Colitis	Colon/rectum infection + circulation support	70	Yes	P&G
18 Paralysis	Nerve Recovery, Paralysis Support	Brainwave range + neuropathy inflammation support	56	Yes	Prl

Mode	Illnesses Targeted	Rationale / Explanation	Total Time in minutes	Auto - Advance	Short Form
19 7-Chakras	Energetic balance, stress, emotional state (non-medical sound healing)	Chakra tone mappings from sound-therapy traditions; intended for relaxation and balance, not a medical frequency set.	42	Yes	Chk
20 Beej Mantra	Subtle energy, mantric resonance (non-medical sound healing).	Seed-syllable (“Beej”) tone mappings; described as experimental/spiritual heuristic practice / treatment.	49	Yes	BJM
21 Brain Wave	Sleep, anxiety, cognition, emotional balance (sound healing).	Uses brainwave ranges (Delta/Theta/Alpha/Beta/Gamma) and Solfeggio-style tones commonly used in binaural/entrainment discussions for calming and focus.	98	Yes	Brn
22 Healing Mode	Pain relief, tissue healing, relaxation, emotional healing (sound healing).	Solfeggio-tradition “healing tones” bundle for mood/relaxation and supportive pain/tissue themes; positioned as non-medical sound healing.	133	Yes	Hea
23 Detox Body	General Detoxification, Cleansing	Solfeggio cleanse + Rife detox core	35	Yes	Det
24 Immune Boost	Immunity, Lymphatic Drainage, Circulation	Lymphatic/immune helper frequencies	35	Yes	Imm
25 Hormonal Care	Thyroid, Adrenal, Reproductive Hormones	Endocrine infection support + speculative hormone tones	70	Yes	Hor
26 Parasite Sweep	Broad Parasite Spectrum (Worms, Flukes, Protozoa)	Targeted parasite range sweep; covers common parasite harmonics	30	Continuous	PSw

Mode	Illnesses Targeted	Rationale / Explanation	Total Time in minutes	Auto - Advance	Short Form
27 Pets & Food	Pet Support, Food/Water Energizing (sound healing)	Broad-spectrum frequency for pet wellness and food/water vibrational treatment UP/DOWN selects duration (default 15 min)	5, 10, 15, 20, 25, 30 min	Fixed	P&F
28 Sleep & Relax	Stress relief, Restful slumber (sound healing)	Progressive brainwave entrainment: Alpha → Delta + deep relaxation	55	Yes	S&R
29 Rife Classics	BX virus, staph/strep, foundational pathogens	Original Rife lab 44 frequencies - 7 min each frequency	308	Yes	Rif
30 Factory Reset	Erases flash memory (menu + frequency)	Clears saved settings (last mode + saved manual frequency) and restarts the device.	N/A	N/A	FRE
31 Manual Mode Hz	Custom Single Frequency in Hz between 2 Hz to 999 Hz	User-defined; with default 50 Hz	7	Fixed	MHz
32 Pemphigus Vul	Autoimmune skin blistering support (Pemphigus vulgaris support)	Condition-focused supportive set intended for autoimmune/skin healing support (complementary).	98	Yes	PeV
33 UTI Kidney	UTI / E. coli urinary support, kidney/urinary discomfort support.	Focused urinary-tract infection support set, meant to simplify selection vs general modes.	108	Yes	UTI
34 Skin Wound	Skin/wound infection support.	Condition-focused infection + skin healing support set aimed at common wound/skin infection patterns.	70	Yes	SKN
35 Anemia Care	Anemia / blood support.	Long supportive set intended for blood/iron/anemia-type support themes (complementary).	129	Yes	ANM
36 Cervix HPV	Cervical HPV support.	Focused cervix/HPV supportive set designed for quicker	112	Yes	HPV

Mode	Illnesses Targeted	Rationale / Explanation	Total Time in minutes	Auto - Advance	Short Form
		selection than general reproductive mode.			
37 Menstrual Care	Menstrual issues, endometriosis/fibroid support.	Reproductive / uterine supportive frequency group tailored to menstrual irregularity and fibroid / endometriosis-type support themes.	126	Yes	MNC
38 Vaginitis	Vaginitis/vaginosis/Candida support.	Focused vaginal flora / yeast - support frequency grouping; positioned as supportive.	154	Yes	VGN
39 PID Adnexitis	PID/adnexitis support.	Short targeted pelvic inflammatory support set for fast selection.	49	Yes	PID
40 Urinary Incont	Urinary control / bladder support.	Supportive bladder/urinary control frequency grouping (symptom / support oriented).	77	Yes	UIC
41 STIs RTIs	STI/RTI support set (sexual / reproductive tract infection support).	Focused reproductive infection support grouping intended to cover common STI / RTI-type support patterns (complementary).	105	Yes	STI
42 Ovarian Cysts	Ovarian cyst support.	Very short targeted ovarian-cyst supportive set for quick use.	21	Yes	OVC
43 Menopause Help	Menopause support.	Symptom/support-oriented menopause frequency grouping (complementary).	112	Yes	MNP
44 Thyroid Care	Thyroid support.	Condition-focused thyroid supportive set (not a replacement for medical thyroid treatment).	77	Yes	THY

Mode	Illnesses Targeted	Rationale / Explanation	Total Time in minutes	Auto - Advance	Short Form
45 Bone Health	Bone regeneration / bone health support.	Supportive bone healing/regeneration frequency set for bone health themes.	112	Yes	BON
46 GallBldr Pain	Gallbladder Pain Support	Condition-focused gallbladder supportive grouping intended for gallbladder discomfort themes.	123	Yes	GAL
47 Migraine Relief	Migraine / Headache relief	Headache support	35	Yes	MIG
48 Sinus / Vertigo	Sinusitis (acute/chronic), rhinitis/runny nose, tinnitus (ringing), vertigo/giddiness, ear pressure/ENT congestion themes.	Focused ENT subset combining sinus + nasal + inner-ear support frequencies commonly repeated across broader sets, intended as a more targeted option than general ENT when symptoms are mainly sinus blockage, nasal irritation, ringing, or dizziness.	186	Yes	SRV
49 Eyes Support	Cataracts, Glaucoma, Macular Degeneration, Eye Strain, Dry Eyes, Diabetic Retinopathy, Conjunctivitis, Blurred vision	Eyes Support is an frequency set intended to provide general eye-comfort support (surface irritation/strain) and broad eye-health support themes (lens/retina/optic nerve)	186	Yes	Eye
50 Liver Support	Hangover, Jaundice, Lipoma, Listerioses, Liver stores/bile/cholesterol support, Liver enlargement, Hepatitis A/B/C (including chronic, cryptogenic, viral), Liver cirrhosis, Liver diseases, Fatty liver, Parasites flukes liver, Pesticide detox, Surgery	Comprehensive liver health and detoxification support covering multiple hepatitis types, cirrhosis, fatty liver, jaundice, hangover recovery, pesticide detoxification, parasitic liver flukes, and post-anesthesia recovery. Complementary to medical treatment, not a replacement. For	116	Yes	Liv

Mode	Illnesses Targeted	Rationale / Explanation	Total Time in minutes	Auto - Advance	Short Form
	anesthesia detox, Vitiligo, Yellow fever	acute liver conditions, seek immediate medical care.			
51 Manual Mode KHz	Custom Single Frequency in KHz between 1 KHz to 500 KHz	User-defined; with default 30 KHz of Clark Hulda univeral frequency.	7	No	MKH

Detailed Sweep & Manual Mode Explanations

Mode 03: Sweep Mode

- **Formula**: 2 Hz → 500 kHz over 60 minutes
- "Full spectrum sweeps" are a **classic technique** for unknown pathogens.
- **Rationale**: Covers **all known pathogen ranges**:

Range	Targets
2-1000 Hz	Pain, Inflammation, Low-Frequency Bacteria
1000-10kHz	Most Viruses, Common Bacteria, Candida
10k-100kHz	Parasites, Flukes, Advanced Infections
100k-500kHz	Experimental High-Range Pathogens

Mode 26: Parasite Sweep

- **Formula**: ~20 Hz → 50 kHz over 30 minutes
- Parasite frequencies cluster between **400 Hz-50 kHz** (tapeworms 15-40kHz, flukes 20-30kHz).
- **Rationale**: More targeted than full sweep; focuses on documented parasite harmonic ranges while avoiding ultra-low/high extremes.

Mode 31: Manual Mode Hz

- **Formula**: Starts at 50 Hz; changes by +25 Hz / −1 Hz below (2 Hz minimum, 999 Hz maximum)
- It allows custom frequency in Hz to be set for 7 min.

Mode 51: Manual Mode KHz

- **Formula**: Starts at 30 KHz (Dr. Clark standard); changes by +25 KHz / −1 KHz below (1 KHz minimum, 500 KHz maximum)
- 30 kHz is the **Hulda Clark zapper frequency** – specifically chosen because many pathogens resonate/die at this square-wave frequency.
- **Rationale**: Allows experimentation with exact frequencies from Rife lists while maintaining 7-minute session safety timer. Perfect for targeting specific pathogens identified in research.
- It allows custom frequency in KHz to be set for 7 min.

Special Modes (Non-Medical / Sound Healing)

Sound Healing Modes (19-22, 27, 28): For Chakras, Brain Waves, Healing, and Pets & Food, use the External Speaker (with volume control) with AUX Cable instead of electrodes.



These gentle frequencies work through sound vibrations. Many are too low to hear but still help relax your mind and body. Connect the External Speaker with AUX Cable, play through External Speaker with AUX Cable, and let the calming tones do their work.

Mode	Focus	Source	Explanation	Total Duration in minutes	Auto-Advance
19 7-Chakras	Energy Center Balance	Sound healing	Chakra tone mappings	42	Yes
20 Beej Mantra	Subtle Energy / Mantras	Spiritual heuristics	Seed syllable frequencies; experimental	49	Yes
21 Brain Wave	Sleep, Anxiety, Cognition	Binaural + Solfeggio	Brainwave entrainment + emotional healing tones	98	Yes
22 Healing Mode	Pain, Tissue Repair, Emotional Healing	Solfeggio tradition	Extended healing frequency set; mood/dopamine support	112	Yes
27 Pets & Food	Pet Support, Food Energizing	General support	Broad-spectrum frequency for pet wellness and food/water vibrational treatment	5 to 30	Fixed
28 Sleep & Relax	Stress relief, Restful slumber	Sound healing	Progressive brainwave entrainment: Alpha → Delta + deep relaxation	55 min	Yes



01 Program Mode

- **Dr. Piyush Saxena's Protocol**: These exact frequencies (2128 Hz, 10 kHz, 30 kHz, 48 kHz) form his core "Cleansing Therapy" sequence for general detoxification, pathogen clearance, and systemic reset. Dr. Saxena promotes this as a foundational 4-step program for natural healing without drugs, often run as 7-minute steps with breaks.

- **Standards**:
  - **2128 Hz**: Ubiquitous "mortal oscillatory rate" (MOR) for broad-spectrum pathogens (bacteria, viruses, abdominal infections).
  - **10 kHz**: General antimicrobial, detoxification support (Dr. Clark/Hulda Clark influence).
  - **30 kHz**: Dr. Hulda Clark's universal parasite/virus killer (7-20-7-20-7 protocol at 30 kHz).
  - **48 kHz**: High-frequency systemic cleanse, cancer support, and deep pathogen targeting.

Protocol Structure: Each runs for 7 minutes (total ~28 minutes), auto-advancing. UP/DOWN skips steps. This matches empirical Rife usage (e.g., 45+ min on 2128 Hz for cancer) but is optimized for daily use.

Scientific Context: Rife therapy targets microbial resonance via electromagnetic fields. While anecdotal/clinical, no mainstream medical validation exists. Cancer Research UK notes no cancer cure evidence. Use as complementary only.

Detailed Illness Associations and Frequency Explanations

Step 1: 2128 Hz (Dr. Piyush Saxena's Primary Frequency)

Targeted Illnesses/Conditions:

- **General Infections**: Abdominal infections, chlamydia, cytomegalovirus (CMV), ear disorders, intestinal disorders.
- **Cancer Support**: Breast cancer (primary/secondary), prostate cancer, sarcoma, melanoma (use 56+ min).
- **Prostatitis/BPH**: 10+ min essential.
- **Chronic Conditions**: Allergies, E. coli, hookworm, hypertension support.

Mechanism: Rife's "tertiary cancer frequency" and broad-spectrum killer. Disrupts pathogen cell walls via resonance. Dr. Saxena starts here for foundational cleansing.

Step 2: 10 kHz (Systemic Antimicrobial)

Targeted Illnesses/Conditions:

- **Athlete's Foot, Aura Cleansing**: Fungal/mold support.
- **Dental Infections, Dandruff**: Secondary oral/skin pathogens.
- **General Detox**: Bacteria, prostatitis, colorectal issues.
- **Cancer Adjunct**: Prostate, general carcinoma.

Mechanism: High-energy frequency for biofilm disruption and deep-tissue penetration. Complements 2128 Hz for multi-layered pathogen elimination.

Step 3: 30 kHz (Dr. Hulda Clark Universal Killer)

Targeted Illnesses/Conditions:

- **Parasites/Viruses**: All parasites (Worms, Flukes, Protozoa), viruses (EBV, HIV support), Bacteria Generally.
- **Liver/Gallbladder Cleanse**: Dr. Saxena's heart attack prevention protocol post-cleanse.
- **Immune Boost**: Typhoid, Klebsiella, General infections.

Mechanism: Hulda Clark's single-frequency protocol (30 kHz for 7 min x3 with 20-min breaks) destroys all Pathogens via Square-Wave Zapping. Simplifies sets; Dr. Saxena adapts for daily use.

Step 4: 48 kHz (Deep Systemic Reset)

Targeted Illnesses/Conditions:

- **Cancer (Advanced)**: BX cancer virus (use with 21275 Hz if capable), lymphoma support.
- **High-Frequency Pathogens**: Resistant strains, deep-tissue infections not hit by lower freqs.
- **Prostate/Colorectal**: Final cleanse phase.

Mechanism: Ultra-high frequency for "finishing" resistant microbes. Targets high-resonance pathogens missed by lower bands. Completes Saxena's protocol for full-spectrum coverage.

Usage Recommendations for Program Mode

- **Daily**: Run full 28-min cycle 1-2x/day for maintenance.
- **Intensive**: Extend 2128 Hz (UP to repeat) for cancer/chronic issues.
- **Integration**: Follow with Common Toxins or Parasite Care modes.
- **Disclaimer**: Experimental/complementary. Not FDA-approved; consult physician. No proven cure for diseases like cancer.

This validation confirms the frequencies align with Dr. Saxena and Rife-style frequency traditions, based on multiple independent references and commonly used frequency lists.

30 Factory Reset

Focus: If you want to reset the Zenthera Zapper to restore Manual Mode back to 30 kHz and clear the last-used menu selection, use the Factory Reset option. This erases saved last mode + last manual frequency and restarts the device.

Steps:

- From the mode selection menu, navigate to “**30 Factory Reset**” and press MODE to enter the reset confirmation screen.
- The device will display: “**Factory Reset**?” and “**UP:Y DOWN/MODE:N**”.
- Press UP button to confirm (triple-beep confirmation) and perform the Reset + Reboot, or press DOWN or MODE button to Cancel and return to the menu.
- Factory reset does not change hardware, only clears memory/settings.

Women’s Health

Note: In addition to broad modes, the firmware includes dedicated women’s health and Condition-focused modes for faster selection:

33 UTI Kidney
Illnesses/Focus: UTI / urinary tract support (E. coli–focused set).

34 Skin Wound
Illnesses/Focus: Skin/wound infection support.

35 Anemia Care
Illnesses/Focus: Anemia / blood support.

36 Cervix HPV
Illnesses/Focus: Cervical HPV support.

37 Menstrual Care
Illnesses/Focus: Menstrual disorders; endometriosis and fibroid/cyst support.

38 Vaginitis
Illnesses/Focus: Vaginitis/vaginosis; Candida/yeast support.

39 PID Adnexitis
Illnesses/Focus: PID (pelvic inflammatory disease) / adnexitis support.

40 Urinary Incont
Illnesses/Focus: Urinary incontinence; bladder/prostate/urinary control support.

41 STIs RTIs
Illnesses/Focus: STI/RTI support (includes gonorrhea/chlamydia/cervicitis-oriented set).

42 Ovarian Cysts
Illnesses/Focus: Ovarian cyst support.

43 Menopause Help
Illnesses/Focus: Menopause support.

Other Specialized Health Modes

44 Thyroid Care
Illnesses/Focus: Thyroid support (hypothyroidism, hyperthyroidism, thyroiditis).

45 Bone Health
Illnesses/Focus: Bone regeneration; fractures/osteoporosis/bone healing support.

46 GallBldr Pain
Illnesses/Focus: Gallbladder support (gallbladder pain/gallstones-oriented set).

47 Migraine Relief
Illnesses/Focus: Migraine / Headache relief.

48 Sinus/Vertigo Focus
Illnesses/Focus: Sinusitis (acute/chronic), rhinitis/runny nose, tinnitus (ringing), vertigo/giddiness, ear pressure/ENT congestion themes.

49 Eyes Support

Illnesses/Focus: Cataracts, Macular Degeneration, Glaucoma, Eye Strain, Dry Eyes, Diabetic Retinopathy, Conjunctivitis, Blurred vision

50 Liver Support

Illnesses/Focus: Comprehensive liver health and detoxification support covering multiple hepatitis types, cirrhosis, fatty liver, jaundice, hangover recovery, pesticide detoxification, parasitic liver flukes, and post-anesthesia recovery. Must run **14 Parasite Care along with 50 Liver Support**.

Suggested Mode Reference

This table is a wellness-support lookup for choosing modes; it is not diagnostic, and results are not guaranteed. For emergencies or serious symptoms, seek medical care.

Note: **Mode 02 Core Support** has list of frequencies to support most of the Virus, Parasite, Cancer, Illness, Bacteria, Pathogen, Condition, Fungus etc.

S/No	Condition / Pathogen	Type	Primary mode	Secondary mode(s)
1	Acanthamoeba infection	Parasite	14 Parasite Care	02 Core Support
2	Acinetobacter baumannii infection	Bacteria	04 Common Toxins	02 Core Support
3	Acne	Condition	34 Skin Wound	04 Common Toxins, 02 Core Support
4	Actinomyces infection	Bacteria	04 Common Toxins	02 Core Support
5	Acute bronchitis	Condition	10 Organ Health	04 Common Toxins, 02 Core Support
6	Acute gastroenteritis	Condition	16 Ulcer Care	04 Common Toxins, 02 Core Support
7	Acute hepatitis	Condition	50 Liver Support	05 Viral Relief, 02 Core Support
8	Acute lymphoblastic leukemia	Cancer	15 Cancer Care	02 Core Support
9	Acute myeloid leukemia	Cancer	15 Cancer Care	02 Core Support
10	Acute pancreatitis	Condition	10 Organ Health	02 Core Support
11	Acute sinusitis	Condition	48 Sinus/Vertigo	06 ENT Care, 04 Common Toxins, 02 Core Support

12	Acute viral hepatitis (unspecified)	Virus	50 Liver Support	05 Viral Relief, 02 Core Support
13	Acute viral respiratory infection	Virus	05 Viral Relief	04 Common Toxins, 02 Core Support
14	Adenovirus conjunctivitis	Virus	49 Eyes Support	05 Viral Relief, 02 Core Support
15	Adenovirus gastroenteritis	Virus	05 Viral Relief	16 Ulcer Care, 02 Core Support
16	Adenovirus respiratory	Virus	05 Viral Relief	02 Core Support
17	Adrenal cancer	Cancer	15 Cancer Care	02 Core Support
18	Aeromonas infection	Bacteria	04 Common Toxins	02 Core Support
19	Allergic rhinitis	Condition	48 Sinus/Vertigo	06 ENT Care, 02 Core Support
20	Amoeba liver abscess	Parasite	14 Parasite Care	50 Liver Support, 02 Core Support
21	Amoebiasis	Parasite	14 Parasite Care	02 Core Support
22	Anemia (iron deficiency)	Condition	35 Anemia Care	02 Core Support
23	Angina pectoris	Condition	08 Cardiac Care	02 Core Support
24	Anisakiasis	Parasite	14 Parasite Care	02 Core Support
25	Anthrax	Bacteria	04 Common Toxins	02 Core Support
26	Anxiety disorder	Condition	22 Healing Mode	21 Brain Wave, 28 Sleep & Relax
27	Aphthous ulcers	Condition	16 Ulcer Care	07 Dental Care, 02 Core Support
28	Arrhythmia	Condition	08 Cardiac Care	02 Core Support
29	Arthritis (osteo)	Condition	13 Arthritis Care	22 Healing Mode, 02 Core Support
30	Ascariasis	Parasite	14 Parasite Care	02 Core Support
31	Aspergillosis	Fungus	04 Common Toxins	02 Core Support
32	Asthma (chronic)	Condition	04 Common Toxins	10 Organ Health, 02 Core Support
33	Atrial fibrillation	Condition	08 Cardiac Care	02 Core Support
34	Avian influenza (H5N1 risk)	Virus	05 Viral Relief	02 Core Support
35	Babesiosis (rare)	Parasite	14 Parasite Care	02 Core Support
36	Bacillus cereus food poisoning	Bacteria	04 Common Toxins	17 Piles & Gut, 02 Core Support
37	Bacterial conjunctivitis	Bacteria	49 Eyes Support	04 Common Toxins, 02 Core Support
38	Bacterial infection (bloodstream, unspecified)	Bacteria	04 Common Toxins	02 Core Support
39	Bacterial infection (gastrointestinal, unspecified)	Bacteria	04 Common Toxins	17 Piles & Gut, 02 Core Support
40	Bacterial infection (post-surgical, unspecified)	Bacteria	04 Common Toxins	34 Skin Wound, 02 Core Support

41	Bacterial infection (respiratory, unspecified)	Bacteria	04 Common Toxins	10 Organ Health, 02 Core Support
42	Bacterial infection (skin, unspecified)	Bacteria	34 Skin Wound	04 Common Toxins, 02 Core Support
43	Bacterial infection (urinary, unspecified)	Bacteria	33 UTI Kidney	04 Common Toxins, 02 Core Support
44	Bacterial infection (wound, unspecified)	Bacteria	34 Skin Wound	02 Core Support
45	Bacterial meningitis	Bacteria	04 Common Toxins	02 Core Support
46	Bacterial vaginosis	Bacteria	38 Vaginitis	41 STIs RTIs, 02 Core Support
47	Bartonella henselae (cat-scratch)	Bacteria	04 Common Toxins	02 Core Support
48	Basal cell carcinoma	Cancer	15 Cancer Care	02 Core Support
49	Bedbug infestation	Parasite	14 Parasite Care	02 Core Support
50	Benign prostatic hyperplasia	Condition	11 Reproductive Organs	25 Hormonal Care, 02 Core Support
51	BK virus	Virus	05 Viral Relief	02 Core Support
52	Bladder cancer	Cancer	15 Cancer Care	02 Core Support
53	Blastomycosis (rare)	Fungus	04 Common Toxins	02 Core Support
54	Blurred vision	Condition	49 Eyes Support	10 Organ Health, 02 Core Support
55	Boils / furunculosis	Condition	34 Skin Wound	04 Common Toxins, 02 Core Support
56	Bronchiectasis	Condition	10 Organ Health	04 Common Toxins, 02 Core Support
57	Brucella melitensis infection	Bacteria	04 Common Toxins	02 Core Support
58	Brucellosis	Bacteria	04 Common Toxins	02 Core Support
59	Campylobacter enteritis	Bacteria	04 Common Toxins	17 Piles & Gut, 02 Core Support
60	Campylobacter jejuni infection	Bacteria	04 Common Toxins	17 Piles & Gut, 02 Core Support
61	Candidiasis	Fungus	04 Common Toxins	38 Vaginitis, 07 Dental Care, 02 Core Support
62	Cataracts	Condition	49 Eyes Support	10 Organ Health, 02 Core Support
63	Cellulitis	Bacteria	34 Skin Wound	04 Common Toxins, 02 Core Support
64	Cervicitis / STI / RTI (general)	Condition	41 STIs RTIs	36 Cervix HPV, 02 Core Support
65	Chandipura virus	Virus	05 Viral Relief	02 Core Support
66	Chikungunya	Virus	05 Viral Relief	02 Core Support
67	Chlamydia	Bacteria	41 STIs RTIs	11 Reproductive Organs, 02 Core Support

68	Chlamydophila pneumoniae infection	Bacteria	04 Common Toxins	10 Organ Health, 02 Core Support
69	Cholangiocarcinoma	Cancer	15 Cancer Care	50 Liver Support, 02 Core Support
70	Cholera	Bacteria	04 Common Toxins	17 Piles & Gut, 02 Core Support
71	Chronic diarrhea	Condition	17 Piles & Gut	16 Ulcer Care, 02 Core Support
72	Chronic fatigue syndrome	Condition	22 Healing Mode	13 Arthritis Care, 02 Core Support
73	Chronic kidney disease	Condition	09 Metabolic Care	33 UTI Kidney, 02 Core Support
74	Chronic lymphocytic leukemia	Cancer	15 Cancer Care	02 Core Support
75	Chronic myeloid leukemia	Cancer	15 Cancer Care	02 Core Support
76	Chronic pancreatitis	Condition	10 Organ Health	02 Core Support
77	Chronic pelvic pain	Condition	37 Menstrual Care	11 Reproductive Organs, 25 Hormonal Care
78	Chronic prostatitis	Condition	11 Reproductive Organs	33 UTI Kidney, 25 Hormonal Care
79	Chronic sinusitis	Condition	48 Sinus/Vertigo	06 ENT Care, 02 Core Support
80	Chronic tonsillitis	Condition	06 ENT Care	04 Common Toxins, 02 Core Support
81	Cirrhosis	Condition	50 Liver Support	02 Core Support
82	Citrobacter freundii infection	Bacteria	04 Common Toxins	02 Core Support
83	Clonorchiasis (rare)	Parasite	14 Parasite Care	50 Liver Support, 02 Core Support
84	Clostridioides difficile infection	Bacteria	04 Common Toxins	17 Piles & Gut, 02 Core Support
85	Clostridium perfringens food poisoning	Bacteria	04 Common Toxins	17 Piles & Gut, 02 Core Support
86	CMV infection	Virus	05 Viral Relief	02 Core Support
87	Coccidioidomycosis (travel risk)	Fungus	04 Common Toxins	02 Core Support
88	Community-acquired pneumonia	Condition	10 Organ Health	04 Common Toxins, 02 Core Support
89	Constipation (functional)	Condition	17 Piles & Gut	02 Core Support
90	COPD	Condition	10 Organ Health	02 Core Support
91	Coronary artery disease	Condition	08 Cardiac Care	02 Core Support
92	Coronavirus 229E	Virus	05 Viral Relief	02 Core Support
93	Coronavirus HKU1	Virus	05 Viral Relief	02 Core Support
94	Coronavirus NL63	Virus	05 Viral Relief	02 Core Support
95	Coronavirus OC43	Virus	05 Viral Relief	02 Core Support
96	COVID-19	Virus	05 Viral Relief	02 Core Support
97	Coxiella burnetii (Q fever)	Bacteria	04 Common Toxins	02 Core Support

98	Coxsackievirus A16	Virus	05 Viral Relief	02 Core Support
99	Coxsackievirus B	Virus	05 Viral Relief	02 Core Support
100	Crimean-Congo hemorrhagic fever	Virus	05 Viral Relief	02 Core Support
101	Cryptococcosis	Fungus	04 Common Toxins	02 Core Support
102	Cryptosporidiosis	Parasite	14 Parasite Care	02 Core Support
103	Cutaneous larva migrans	Parasite	14 Parasite Care	02 Core Support
104	Cutaneous leishmaniasis	Parasite	14 Parasite Care	02 Core Support
105	Cyclospora infection	Parasite	14 Parasite Care	02 Core Support
106	Cystoisospora infection	Parasite	14 Parasite Care	02 Core Support
107	Dengue	Virus	05 Viral Relief	02 Core Support
108	Dental abscess	Bacteria	07 Dental Care	04 Common Toxins, 02 Core Support
109	Dental caries	Condition	07 Dental Care	02 Core Support
110	Dental fistula (Fistula dentalis)	Condition	16 Ulcer Care	07 Dental Care, 02 Core Support
111	Depression	Condition	22 Healing Mode	21 Brain Wave, 28 Sleep & Relax
112	Dermatitis / eczema	Condition	10 Organ Health	34 Skin Wound, 02 Core Support
113	Dermatophytosis (ringworm)	Fungus	04 Common Toxins	34 Skin Wound, 02 Core Support
114	Diabetic neuropathy	Condition	13 Arthritis Care	22 Healing Mode, 09 Metabolic Care
115	Diabetic retinopathy	Condition	49 Eyes Support	09 Metabolic Care, 10 Organ Health, 02 Core Support
116	Diphtheria	Bacteria	04 Common Toxins	02 Core Support
117	Dirofilariasis (rare)	Parasite	14 Parasite Care	02 Core Support
118	Dracunculiasis	Parasite	14 Parasite Care	02 Core Support
119	Dry eyes	Condition	49 Eyes Support	10 Organ Health, 02 Core Support
120	Dysmenorrhea	Condition	37 Menstrual Care	25 Hormonal Care, 02 Core Support
121	E. coli diarrhea	Bacteria	04 Common Toxins	17 Piles & Gut, 02 Core Support
122	EBV infection	Virus	05 Viral Relief	02 Core Support
123	Echinococcosis (hydatid)	Parasite	14 Parasite Care	02 Core Support
124	Echovirus	Virus	05 Viral Relief	02 Core Support
125	Endometrial cancer	Cancer	15 Cancer Care	02 Core Support
126	Endometriosis	Condition	37 Menstrual Care	11 Reproductive Organs, 25 Hormonal Care
127	Enterobacter cloacae infection	Bacteria	04 Common Toxins	02 Core Support
128	Enterobiasis (pinworm)	Parasite	14 Parasite Care	02 Core Support

129	Enterococcus faecalis infection	Bacteria	33 UTI Kidney	04 Common Toxins, 02 Core Support
130	Enterovirus (other)	Virus	05 Viral Relief	02 Core Support
131	Enterovirus 71	Virus	05 Viral Relief	02 Core Support
132	Enterovirus D68	Virus	05 Viral Relief	02 Core Support
133	Enterovirus infection	Virus	05 Viral Relief	02 Core Support
134	Epilepsy	Condition	21 Brain Wave	28 Sleep & Relax, 02 Core Support
135	Ewing sarcoma	Cancer	15 Cancer Care	02 Core Support
136	Eye strain / computer vision syndrome	Condition	49 Eyes Support	22 Healing Mode, 10 Organ Health
137	Fascioliasis (rare)	Parasite	14 Parasite Care	50 Liver Support, 02 Core Support
138	Fatty liver disease	Condition	50 Liver Support	23 Detox Body, 02 Core Support
139	Fibromyalgia	Condition	13 Arthritis Care	22 Healing Mode, 02 Core Support
140	Filariasis hydrocele	Parasite	14 Parasite Care	02 Core Support
141	Fistula (general)	Condition	16 Ulcer Care	02 Core Support
142	Flea infestation	Parasite	14 Parasite Care	02 Core Support
143	Food poisoning (bacterial)	Bacteria	04 Common Toxins	17 Piles & Gut, 02 Core Support
144	Francisella tularensis (tularemia risk)	Bacteria	04 Common Toxins	02 Core Support
145	Frozen shoulder	Condition	13 Arthritis Care	22 Healing Mode, 02 Core Support
146	Fungal otitis externa	Fungus	06 ENT Care	04 Common Toxins, 02 Core Support
147	Fungal sinusitis	Fungus	48 Sinus/Vertigo	06 ENT Care, 02 Core Support
148	Fungal skin infection (unspecified)	Fungus	34 Skin Wound	04 Common Toxins, 02 Core Support
149	Fungal vaginitis	Fungus	38 Vaginitis	02 Core Support
150	Gallbladder cancer	Cancer	15 Cancer Care	46 GallBldr Pain, 02 Core Support
151	Gallbladder pain / gallstones support	Condition	46 GallBldr Pain	50 Liver Support, 02 Core Support
152	Gastritis (bacterial)	Bacteria	16 Ulcer Care	04 Common Toxins, 02 Core Support
153	GERD	Condition	10 Organ Health	16 Ulcer Care, 02 Core Support
154	Giardia chronic	Parasite	14 Parasite Care	02 Core Support
155	Giardiasis	Parasite	14 Parasite Care	02 Core Support
156	Gingivitis	Condition	07 Dental Care	02 Core Support
157	GIST (gastrointestinal stromal tumor)	Cancer	15 Cancer Care	02 Core Support
158	Glaucoma	Condition	49 Eyes Support	10 Organ Health, 02 Core Support
159	Glioblastoma	Cancer	15 Cancer Care	02 Core Support
160	Glioma (other)	Cancer	15 Cancer Care	02 Core Support

161	Gnathostomiasis	Parasite	14 Parasite Care	02 Core Support
162	Gonorrhea	Bacteria	41 STIs RTIs	11 Reproductive Organs, 02 Core Support
163	Gout	Condition	13 Arthritis Care	02 Core Support
164	Graves disease	Condition	44 Thyroid Care	25 Hormonal Care, 02 Core Support
165	Group B strep infection	Bacteria	04 Common Toxins	02 Core Support
166	H. pylori infection	Bacteria	16 Ulcer Care	04 Common Toxins, 02 Core Support
167	Haemophilus ducreyi (chancroid)	Bacteria	41 STIs RTIs	02 Core Support
168	Haemophilus influenzae infection	Bacteria	04 Common Toxins	10 Organ Health, 02 Core Support
169	Haemorrhoids	Condition	17 Piles & Gut	02 Core Support
170	Hand foot mouth disease	Virus	05 Viral Relief	02 Core Support
171	Hangover	Condition	50 Liver Support	23 Detox Body
172	Hantavirus infection	Virus	05 Viral Relief	02 Core Support
173	Hashimoto thyroiditis	Condition	44 Thyroid Care	12 Autoimmune, 25 Hormonal Care, 02 Core Support
174	Headache (tension / cluster / sinus)	Condition	47 Migraine Relief	48 Sinus/Vertigo, 22 Healing Mode, 02 Core Support
175	Heart failure	Condition	08 Cardiac Care	02 Core Support
176	Helicobacter-related ulcer	Bacteria	16 Ulcer Care	04 Common Toxins, 02 Core Support
177	Helminth infection (unspecified)	Parasite	14 Parasite Care	02 Core Support
178	Hepatitis A	Virus	50 Liver Support	05 Viral Relief, 02 Core Support
179	Hepatitis B	Virus	50 Liver Support	05 Viral Relief, 02 Core Support
180	Hepatitis C	Virus	50 Liver Support	05 Viral Relief, 02 Core Support
181	Hepatitis D	Virus	50 Liver Support	05 Viral Relief, 02 Core Support
182	Hepatitis E	Virus	50 Liver Support	05 Viral Relief, 02 Core Support
183	Hepatitis G / GBV-C	Virus	50 Liver Support	05 Viral Relief, 02 Core Support
184	Herpes simplex	Virus	05 Viral Relief	02 Core Support
185	Herpes zoster	Virus	05 Viral Relief	22 Healing Mode, 02 Core Support
186	HHV-6	Virus	05 Viral Relief	02 Core Support
187	HHV-7	Virus	05 Viral Relief	02 Core Support
188	HHV-8 (KSHV)	Virus	05 Viral Relief	02 Core Support
189	Histoplasmosis	Fungus	04 Common Toxins	02 Core Support
190	HIV	Virus	05 Viral Relief	02 Core Support

191	Hodgkin lymphoma	Cancer	15 Cancer Care	02 Core Support
192	Hookworm disease	Parasite	14 Parasite Care	02 Core Support
193	Hospital-acquired pneumonia	Condition	10 Organ Health	04 Common Toxins, 02 Core Support
194	HPV infection	Virus	36 Cervix HPV	41 STIs RTIs, 11 Reproductive Organs
195	HTLV-1	Virus	05 Viral Relief	02 Core Support
196	Human astrovirus	Virus	05 Viral Relief	02 Core Support
197	Human bocavirus	Virus	05 Viral Relief	02 Core Support
198	Human metapneumovirus	Virus	05 Viral Relief	02 Core Support
199	Human parechovirus	Virus	05 Viral Relief	02 Core Support
200	Hydatid disease	Parasite	14 Parasite Care	02 Core Support
201	Hymenolepiasis	Parasite	14 Parasite Care	02 Core Support
202	Hypertension	Condition	08 Cardiac Care	02 Core Support
203	Hypertensive urgency	Condition	08 Cardiac Care	02 Core Support
204	Hyperthyroidism	Condition	44 Thyroid Care	25 Hormonal Care, 02 Core Support
205	Hypothyroidism	Condition	44 Thyroid Care	25 Hormonal Care, 02 Core Support
206	IBS	Condition	17 Piles & Gut	16 Ulcer Care, 02 Core Support
207	Influenza A	Virus	05 Viral Relief	02 Core Support
208	Influenza B	Virus	05 Viral Relief	02 Core Support
209	Insomnia	Condition	28 Sleep & Relax	21 Brain Wave, 22 Healing Mode
210	Interstitial lung disease	Condition	10 Organ Health	02 Core Support
211	Intestinal parasites (mixed)	Parasite	14 Parasite Care	02 Core Support
212	Intestinal worm infection (unspecified)	Parasite	14 Parasite Care	02 Core Support
213	Japanese encephalitis	Virus	05 Viral Relief	02 Core Support
214	Jaundice (non-obstructive)	Condition	50 Liver Support	02 Core Support
215	JC virus	Virus	05 Viral Relief	02 Core Support
216	Kala-azar (visceral leishmaniasis)	Parasite	14 Parasite Care	02 Core Support
217	Kidney stones	Condition	09 Metabolic Care	33 UTI Kidney, 02 Core Support
218	Klebsiella granulomatis (donovanosis)	Bacteria	41 STIs RTIs	02 Core Support
219	Klebsiella pneumoniae infection	Bacteria	04 Common Toxins	10 Organ Health, 02 Core Support
220	Kyasanur Forest Disease	Virus	05 Viral Relief	02 Core Support
221	Laryngeal cancer	Cancer	15 Cancer Care	02 Core Support
222	Lassa fever (imported risk)	Virus	05 Viral Relief	02 Core Support
223	Legionella infection	Bacteria	04 Common Toxins	10 Organ Health, 02 Core Support

224	Legionellosis	Bacteria	04 Common Toxins	10 Organ Health, 02 Core Support
225	Leishmaniasis PKDL	Parasite	14 Parasite Care	02 Core Support
226	Leprosy	Bacteria	04 Common Toxins	02 Core Support
227	Leptospira interrogans infection	Bacteria	04 Common Toxins	02 Core Support
228	Leptospirosis	Bacteria	04 Common Toxins	02 Core Support
229	Lipoma	Condition	50 Liver Support	23 Detox Body, 02 Core Support
230	Listeria monocytogenes infection	Bacteria	50 Liver Support	04 Common Toxins, 02 Core Support
231	Liver cirrhosis	Condition	50 Liver Support	02 Core Support
232	Liver enlargement (hepatomegaly)	Condition	50 Liver Support	09 Metabolic Care, 02 Core Support
233	Low back pain	Condition	13 Arthritis Care	22 Healing Mode, 02 Core Support
234	Lyme disease (rare/imported)	Bacteria	14 Parasite Care	04 Common Toxins, 02 Core Support
235	Lymphatic filariasis	Parasite	14 Parasite Care	02 Core Support
236	Macular degeneration (AMD)	Condition	49 Eyes Support	10 Organ Health, 02 Core Support
237	Malaria (P. falciparum)	Parasite	14 Parasite Care	05 Viral Relief, 02 Core Support
238	Malaria (P. vivax)	Parasite	14 Parasite Care	05 Viral Relief, 02 Core Support
239	Malaria relapse / fever (vivax)	Parasite	14 Parasite Care	05 Viral Relief, 02 Core Support
240	Measles	Virus	05 Viral Relief	02 Core Support
241	Melanoma	Cancer	15 Cancer Care	02 Core Support
242	Melioidosis	Bacteria	04 Common Toxins	02 Core Support
243	Meningioma	Cancer	15 Cancer Care	02 Core Support
244	Menopause symptoms	Condition	43 Menopause Help	25 Hormonal Care, 02 Core Support
245	MERS-CoV (imported risk)	Virus	05 Viral Relief	02 Core Support
246	Metabolic syndrome	Condition	09 Metabolic Care	02 Core Support
247	Microsporidiosis	Parasite	14 Parasite Care	02 Core Support
248	Migraine	Condition	47 Migraine Relief	22 Healing Mode, 02 Core Support
249	Morganella morganii infection	Bacteria	04 Common Toxins	02 Core Support
250	Mpox	Virus	05 Viral Relief	02 Core Support
251	MRSA infection	Bacteria	34 Skin Wound	04 Common Toxins, 02 Core Support
252	MSSA infection	Bacteria	34 Skin Wound	04 Common Toxins, 02 Core Support
253	Mucormycosis	Fungus	04 Common Toxins	02 Core Support
254	Multidrug-resistant TB	Bacteria	04 Common Toxins	02 Core Support
255	Multiple myeloma	Cancer	15 Cancer Care	02 Core Support
256	Mumps	Virus	05 Viral Relief	02 Core Support

257	Mycoplasma genitalium infection	Bacteria	41 STIs RTIs	02 Core Support
258	Mycoplasma pneumoniae infection	Bacteria	04 Common Toxins	10 Organ Health, 02 Core Support
259	Myiasis	Parasite	14 Parasite Care	02 Core Support
260	Naegleria fowleri infection	Parasite	14 Parasite Care	02 Core Support
261	Nasopharyngeal carcinoma	Cancer	15 Cancer Care	02 Core Support
262	Neisseria meningitidis infection	Bacteria	04 Common Toxins	02 Core Support
263	Neuroblastoma	Cancer	15 Cancer Care	02 Core Support
264	Neurocysticercosis	Parasite	14 Parasite Care	02 Core Support
265	Nipah	Virus	05 Viral Relief	02 Core Support
266	Nocardia infection	Bacteria	04 Common Toxins	02 Core Support
267	Non-Hodgkin lymphoma	Cancer	15 Cancer Care	02 Core Support
268	Non-typhoidal Salmonella infection	Bacteria	04 Common Toxins	17 Piles & Gut, 02 Core Support
269	Norovirus gastroenteritis	Virus	05 Viral Relief	16 Ulcer Care, 02 Core Support
270	Obesity	Condition	09 Metabolic Care	02 Core Support
271	Onychomycosis	Fungus	04 Common Toxins	02 Core Support
272	Opisthorchiasis (rare)	Parasite	14 Parasite Care	50 Liver Support, 02 Core Support
273	Oral candidiasis	Condition	07 Dental Care	04 Common Toxins, 02 Core Support
274	Oral thrush	Fungus	07 Dental Care	04 Common Toxins, 02 Core Support
275	Orientia tsutsugamushi infection	Bacteria	04 Common Toxins	02 Core Support
276	Oropharyngeal cancer	Cancer	15 Cancer Care	02 Core Support
277	Osteoporosis	Condition	45 Bone Health	02 Core Support
278	Osteosarcoma	Cancer	15 Cancer Care	02 Core Support
279	Otitis externa	Condition	06 ENT Care	02 Core Support
280	Otitis media (bacterial)	Bacteria	06 ENT Care	04 Common Toxins, 02 Core Support
281	Otitis media (recurrent)	Condition	06 ENT Care	02 Core Support
282	Ovarian cysts	Condition	42 Ovarian Cysts	11 Reproductive Organs, 02 Core Support
283	Pancreatic cancer	Cancer	15 Cancer Care	02 Core Support
284	Parainfluenza	Virus	05 Viral Relief	02 Core Support
285	Parainfluenza type 1	Virus	05 Viral Relief	02 Core Support
286	Parainfluenza type 2	Virus	05 Viral Relief	02 Core Support
287	Parainfluenza type 3	Virus	05 Viral Relief	02 Core Support
288	Parainfluenza type 4	Virus	05 Viral Relief	02 Core Support

289	Parasites flukes liver	Parasite	50 Liver Support	14 Parasite Care, 02 Core Support
290	Paratyphoid fever	Bacteria	04 Common Toxins	02 Core Support
291	Parvovirus B19	Virus	05 Viral Relief	02 Core Support
292	PCOS	Condition	11 Reproductive Organs	25 Hormonal Care, 02 Core Support
293	Pediculosis (lice)	Parasite	14 Parasite Care	02 Core Support
294	Pemphigus vulgaris	Condition	32 Pemphigus Vul	12 Autoimmune, 02 Core Support
295	Penile cancer	Cancer	15 Cancer Care	02 Core Support
296	Peptic ulcer disease	Condition	16 Ulcer Care	02 Core Support
297	Periodontal disease	Condition	07 Dental Care	02 Core Support
298	Periodontitis	Bacteria	07 Dental Care	04 Common Toxins, 02 Core Support
299	Peripheral artery disease	Condition	08 Cardiac Care	02 Core Support
300	Pertussis	Bacteria	04 Common Toxins	10 Organ Health, 02 Core Support
301	Pesticide detox	Condition	50 Liver Support	23 Detox Body, 02 Core Support
302	Pharyngitis	Bacteria	06 ENT Care	04 Common Toxins, 02 Core Support
303	Pharyngitis (recurrent)	Condition	06 ENT Care	02 Core Support
304	PID (pelvic inflammatory disease)	Condition	39 PID Adnexitis	41 STIs RTIs, 11 Reproductive Organs
305	Pimples	Condition	34 Skin Wound	04 Common Toxins, 02 Core Support
306	Pituitary tumor	Cancer	15 Cancer Care	02 Core Support
307	Plague	Bacteria	04 Common Toxins	02 Core Support
308	Pneumococcal pneumonia	Bacteria	04 Common Toxins	10 Organ Health, 02 Core Support
309	Pneumocystis jirovecii pneumonia	Fungus	04 Common Toxins	10 Organ Health, 02 Core Support
310	Polio (wild / VDPV)	Virus	05 Viral Relief	02 Core Support
311	Prostatitis (chronic)	Condition	11 Reproductive Organs	33 UTI Kidney, 25 Hormonal Care
312	Proteus mirabilis infection	Bacteria	33 UTI Kidney	04 Common Toxins, 02 Core Support
313	Protozoal diarrhea (unspecified)	Parasite	14 Parasite Care	17 Piles & Gut, 02 Core Support
314	Protozoal infection (unspecified)	Parasite	14 Parasite Care	02 Core Support
315	Providencia infection	Bacteria	04 Common Toxins	02 Core Support
316	Pseudomonas aeruginosa infection	Bacteria	34 Skin Wound	04 Common Toxins, 02 Core Support
317	Pulmonary fibrosis	Condition	10 Organ Health	22 Healing Mode, 02 Core Support
318	Pyelonephritis	Bacteria	33 UTI Kidney	02 Core Support
319	Rabies	Virus	05 Viral Relief	02 Core Support

320	Recurrent UTI	Condition	33 UTI Kidney	11 Reproductive Organs, 02 Core Support
321	Renal cell carcinoma	Cancer	15 Cancer Care	02 Core Support
322	Retinoblastoma	Cancer	15 Cancer Care	49 Eyes Support, 02 Core Support
323	Rheumatic fever	Bacteria	04 Common Toxins	08 Cardiac Care, 02 Core Support
324	Rheumatoid arthritis	Condition	13 Arthritis Care	12 Autoimmune, 02 Core Support
325	Rhinovirus (other types)	Virus	05 Viral Relief	04 Common Toxins, 02 Core Support
326	Rhinovirus common cold	Virus	05 Viral Relief	04 Common Toxins, 02 Core Support
327	Rickettsia conorii infection	Bacteria	04 Common Toxins	02 Core Support
328	Rickettsial fever (other)	Bacteria	04 Common Toxins	02 Core Support
329	Rotavirus B	Virus	05 Viral Relief	16 Ulcer Care, 02 Core Support
330	Rotavirus C	Virus	05 Viral Relief	16 Ulcer Care, 02 Core Support
331	Rotavirus gastroenteritis	Virus	05 Viral Relief	16 Ulcer Care, 02 Core Support
332	RSV	Virus	05 Viral Relief	02 Core Support
333	RSV A	Virus	05 Viral Relief	02 Core Support
334	RSV B	Virus	05 Viral Relief	02 Core Support
335	Rubella	Virus	05 Viral Relief	02 Core Support
336	Sapovirus	Virus	05 Viral Relief	16 Ulcer Care, 02 Core Support
337	SARS (risk)	Virus	05 Viral Relief	02 Core Support
338	Scabies	Parasite	14 Parasite Care	34 Skin Wound, 02 Core Support
339	Scarlet fever	Bacteria	04 Common Toxins	02 Core Support
340	Schistosomiasis (rare)	Parasite	14 Parasite Care	02 Core Support
341	Sciatica	Condition	13 Arthritis Care	22 Healing Mode, 02 Core Support
342	Scrub typhus	Bacteria	04 Common Toxins	02 Core Support
343	Seasonal coronavirus (unspecified)	Virus	05 Viral Relief	02 Core Support
344	Sepsis	Bacteria	04 Common Toxins	02 Core Support
345	Sepsis syndrome	Condition	04 Common Toxins	02 Core Support
346	Serratia marcescens infection	Bacteria	04 Common Toxins	02 Core Support
347	Shigella flexneri infection	Bacteria	04 Common Toxins	17 Piles & Gut, 02 Core Support
348	Shigella sonnei infection	Bacteria	04 Common Toxins	17 Piles & Gut, 02 Core Support
349	Shigellosis	Bacteria	04 Common Toxins	17 Piles & Gut, 02 Core Support
350	Sinusitis (bacterial)	Bacteria	48 Sinus/Vertigo	06 ENT Care, 04 Common Toxins, 02 Core Support

351	Soft tissue sarcoma	Cancer	15 Cancer Care	02 Core Support
352	Soil-transmitted helminth infection	Parasite	14 Parasite Care	02 Core Support
353	Squamous cell skin cancer	Cancer	15 Cancer Care	02 Core Support
354	Staph skin infection	Bacteria	34 Skin Wound	04 Common Toxins, 02 Core Support
355	Stenotrophomonas maltophilia infection	Bacteria	04 Common Toxins	02 Core Support
356	Strep throat	Bacteria	06 ENT Care	04 Common Toxins, 02 Core Support
357	Stroke (hemorrhagic)	Condition	18 Paralysis	13 Arthritis Care, 02 Core Support
358	Stroke (ischemic)	Condition	18 Paralysis	13 Arthritis Care, 02 Core Support
359	Strongyloidiasis	Parasite	14 Parasite Care	02 Core Support
360	Syphilis	Bacteria	41 STIs RTIs	02 Core Support
361	Taenia solium cysticercosis	Parasite	14 Parasite Care	02 Core Support
362	Taeniasis	Parasite	14 Parasite Care	02 Core Support
363	Testicular cancer	Cancer	15 Cancer Care	02 Core Support
364	Tetanus	Bacteria	04 Common Toxins	02 Core Support
365	Thymoma	Cancer	15 Cancer Care	02 Core Support
366	Thyroid cancer	Cancer	15 Cancer Care	44 Thyroid Care, 02 Core Support
367	Tick infestation	Parasite	14 Parasite Care	02 Core Support
368	Tick-borne encephalitis (travel risk)	Virus	05 Viral Relief	02 Core Support
369	Tinea capitis	Fungus	04 Common Toxins	02 Core Support
370	Tinea corporis	Fungus	04 Common Toxins	02 Core Support
371	Tinea cruris	Fungus	04 Common Toxins	02 Core Support
372	Tinea pedis (athlete foot)	Fungus	04 Common Toxins	02 Core Support
373	Tinea versicolor	Fungus	04 Common Toxins	02 Core Support
374	Tonsillitis	Bacteria	06 ENT Care	04 Common Toxins, 02 Core Support
375	Toxocariasis	Parasite	14 Parasite Care	02 Core Support
376	Toxoplasmosis	Parasite	14 Parasite Care	02 Core Support
377	Toxoplasmosis (congenital)	Parasite	14 Parasite Care	02 Core Support
378	Tracheoesophageal fistula	Condition	16 Ulcer Care	02 Core Support
379	Trachoma	Bacteria	49 Eyes Support	04 Common Toxins, 02 Core Support
380	Traveler diarrhea	Condition	17 Piles & Gut	04 Common Toxins, 02 Core Support
381	Treponema pallidum (syphilis)	Bacteria	41 STIs RTIs	02 Core Support
382	Trichinellosis (rare)	Parasite	14 Parasite Care	02 Core Support

383	Trichomoniasis	Parasite	38 Vaginitis	41 STIs RTIs, 02 Core Support
384	Trichuriasis	Parasite	14 Parasite Care	02 Core Support
385	Tuberculosis	Bacteria	04 Common Toxins	02 Core Support
386	Type 1 diabetes	Condition	09 Metabolic Care	02 Core Support
387	Type 2 diabetes	Condition	09 Metabolic Care	02 Core Support
388	Typhoid fever	Bacteria	04 Common Toxins	05 Viral Relief, 02 Core Support
389	Ureaplasma urealyticum infection	Bacteria	41 STIs RTIs	02 Core Support
390	Urinary incontinence / enuresis / bedwetting	Condition	40 Urinary Incont	33 UTI Kidney, 02 Core Support
391	Uterine fibroids	Condition	37 Menstrual Care	11 Reproductive Organs, 25 Hormonal Care
392	UTI (E. coli)	Bacteria	33 UTI Kidney	11 Reproductive Organs, 02 Core Support
393	Vaginal cancer	Cancer	15 Cancer Care	02 Core Support
394	Vaginitis / vaginosis	Condition	38 Vaginitis	41 STIs RTIs, 02 Core Support
395	Varicella (chickenpox)	Virus	05 Viral Relief	02 Core Support
396	Varicose veins	Condition	22 Healing Mode	08 Cardiac Care, 02 Core Support
397	Venous insufficiency	Condition	08 Cardiac Care	02 Core Support
398	Ventilator-associated pneumonia	Condition	10 Organ Health	04 Common Toxins, 02 Core Support
399	Vibrio parahaemolyticus infection	Bacteria	04 Common Toxins	17 Piles & Gut, 02 Core Support
400	Vibrio vulnificus infection	Bacteria	04 Common Toxins	34 Skin Wound, 02 Core Support
401	Viral bronchiolitis	Virus	05 Viral Relief	10 Organ Health, 02 Core Support
402	Viral conjunctivitis	Virus	49 Eyes Support	05 Viral Relief, 02 Core Support
403	Viral encephalitis (unspecified)	Virus	05 Viral Relief	02 Core Support
404	Viral exanthem (unspecified)	Virus	05 Viral Relief	02 Core Support
405	Viral fever (unspecified)	Virus	05 Viral Relief	04 Common Toxins, 02 Core Support
406	Viral gastroenteritis (unspecified)	Virus	05 Viral Relief	16 Ulcer Care, 02 Core Support
407	Viral meningitis (enteroviral)	Virus	05 Viral Relief	02 Core Support
408	Viral myocarditis (enteroviral)	Virus	05 Viral Relief	08 Cardiac Care, 02 Core Support
409	Vitamin D deficiency	Condition	02 Core Support	23 Detox Body
410	Vitiligo	Condition	12 Autoimmune	50 Liver Support, 02 Core Support
411	VRE (Enterococcus faecium)	Bacteria	04 Common Toxins	02 Core Support
412	Vulvar cancer	Cancer	15 Cancer Care	02 Core Support

413	West Nile fever	Virus	05 Viral Relief	02 Core Support
414	Wilms tumor	Cancer	15 Cancer Care	02 Core Support
415	Yellow fever	Virus	50 Liver Support	05 Viral Relief, 02 Core Support
416	Yersinia enterocolitica infection	Bacteria	04 Common Toxins	17 Piles & Gut, 02 Core Support
417	Zika	Virus	05 Viral Relief	02 Core Support

Product Warranty (Hardware)

Limited Hardware Warranty

The manufacturer warrants that the Zenthera Zapper hardware (circuit board, LCD, buttons, enclosure, and buzzer) will be free from defects in materials and workmanship under normal use for a period of **one (1) year** from the date of purchase.

The power adapter, cables, and electrodes are covered by a limited warranty for **three (3) months** from the date of purchase.

Conductive gel is a consumable item and is provided **without warranty**. Once consumed, it can be purchased separately from nearby medical supply shops.

This warranty covers:

- Manufacturing defects in electronic components
- LCD display failures not caused by physical damage
- Button malfunction due to manufacturing defects
- Power circuit failures under normal operating conditions

This warranty does NOT cover:

- Damage from misuse, abuse, or improper handling
- Water damage or exposure to moisture beyond specified IP rating
- Battery deterioration (batteries are consumable items)
- Normal wear and tear
- Damage from unauthorized modifications or repairs
- Use with incorrect voltage or current specifications

Warranty Claim Process:

Contact the volunteer or our website www.thetempleofhealing.org or write a email on CureYourSelfIndia@gmail.com with description of defect. We will coordinate with you and will repair or replace defective units at its sole discretion.

Disclaimer of Medical Warranty

NO MEDICAL WARRANTY

This device is provided "AS IS" for experimental and research purposes only. The manufacturer makes **NO WARRANTY, EXPRESS OR IMPLIED**, regarding:

- The device's effectiveness for any health condition
- The safety or efficacy of frequency-based therapy
- Medical outcomes or health improvements
- Accuracy of frequency programs for any stated purpose

The manufacturer explicitly disclaims:

- Any warranty of merchantability for medical use
- Any warranty of fitness for a particular medical purpose
- Any claims that this device can diagnose, treat, cure, or prevent disease

This device is **NOT** a medical device. It has **NOT** been evaluated or approved by any regulatory health authority (DCGI, CDSCO, NRA, FDA, CE, TGA, or equivalent).

Limitation of Liability

LIMITATION OF LIABILITY

To the maximum extent permitted by law, the manufacturer's liability is limited to:

- The purchase price of the device
- Repair or replacement of defective hardware (within warranty period)

The manufacturer shall **NOT** be liable for:

- Personal injury or health complications arising from device use
- Consequential, incidental, special, or punitive damages
- Medical expenses or healthcare costs
- Loss of time, inconvenience, or loss of use
- Reliance on frequency programs for health decisions
- Interference with medical devices (including pacemakers)

USER ASSUMES ALL RISK associated with using this experimental device.

Safety Disclaimer (Critical)

SAFETY WARNING & DISCLAIMER

🚨 **DO NOT USE if you have**:

- A pacemaker or any implanted electronic medical device
- Are pregnant or suspect pregnancy

- Have epilepsy or seizure disorders
- Have a heart condition without physician clearance

This device generates electromagnetic frequencies. The long-term health effects are unknown. Use is entirely at your own risk.

The manufacturer provides this device for educational and experimental purposes. Users are solely responsible for their health decisions and should consult qualified healthcare professionals before using frequency-based therapies.

Software/Firmware Disclaimer

FIRMWARE & FREQUENCY PROGRAMS DISCLAIMER

The frequency programs stored in this device are compiled from various public-domain sources and historical Rife frequency lists. The manufacturer makes **NO REPRESENTATION** regarding:

- The accuracy of frequency values
- The safety of any specific frequency
- The therapeutic value of any program
- The scientific validity of frequency-based approaches

Frequency programs may be updated via firmware. Users accept all risks associated with firmware updates and frequency program changes.

No Endorsement Disclaimer

NO MEDICAL ENDORSEMENT

References to health conditions, symptoms, or organs in mode labels (e.g., "Cancer Support," "Cardiac Care," "UTI Kidney") are:

- Descriptive labels for frequency list organization only
- NOT medical claims or treatment promises
- NOT endorsed by medical professionals
- NOT validated by clinical studies

These labels reflect historical Rife frequency list categories and do NOT constitute medical advice or treatment recommendations.

<!- End of File -->