

Zenthera Zapper v1.2 User Guide

Dated : 05 Jun 2026

Table of Contents


| | |
|---|----|
| Discover the Power of Zenthera Zapper | 5 |
| Safety Guidelines..... | 6 |
| Battery Care Instructions | 6 |
| Airline Travel Precautions | 7 |
| Battery Status Indicators | 7 |
| Charging Instructions..... | 7 |
| Pro Tips..... | 8 |
| User Operating Manual..... | 8 |
| Welcome to Zenthera Zapper | 8 |
| Quick Start Guide | 9 |
| 1. Power On..... | 9 |
| 2. Select Mode..... | 9 |
| 3. Start Session | 9 |
| 4. During Session..... | 9 |
| Display Explained..... | 10 |
| Using Carbon Electrodes (Pads / Straps)..... | 10 |
| Core Modes Summary | 11 |
| Daily Essentials (Use Electrodes) | 11 |
| Specialized Health (Use Electrodes)..... | 12 |
| Additional Health Modes (Use Electrodes) | 12 |
| Sound Healing Modes (Use External Speaker with AUX Cable - No Electrodes) | 13 |
| Special Modes..... | 14 |

| | |
|---|----|
| 03 Sweep Mode (60 min) | 14 |
| 26 ParasiteSweep (30 min)..... | 14 |
| 31 Manual Mode (7 min) | 14 |
| Button Functions by State | 14 |
| Complete Session Flow Example..... | 15 |
| Daily Usage Recommendations..... | 15 |
| Beginner (30 min/day)..... | 15 |
| Regular / Maintenance (45-120 min/day) | 15 |
| Intensive Protocol..... | 15 |
| External Speaker with AUX Cable and Volume Control Usage (Sound Healing Modes 19-22, 27, 28)..... | 16 |
| Selecting the external Speakers for Sound Healing Modes..... | 16 |
| Must-have specifications | 17 |
| Recommended (best results) specs | 17 |
| What to avoid..... | 17 |
| Troubleshooting | 17 |
| Frequency Table by Mode | 18 |
| Detailed Sweep & Manual Mode Explanations | 24 |
| Mode 03: Sweep Mode | 24 |
| Mode 26: Parasite Sweep..... | 24 |
| Mode 31: Manual Mode | 24 |
| Special Modes (Non-Medical / Sound Healing) | 25 |
| 01 Program Mode..... | 26 |
| Detailed Illness Associations and Frequency Explanations | 27 |
| Step 1: 2128 Hz (Dr. Piyush Saxena's Primary Frequency) | 27 |
| Step 2: 10 kHz (Systemic Antimicrobial) | 27 |
| Step 3: 30 kHz (Dr. Hulda Clark Universal Killer) | 28 |
| Step 4: 48 kHz (Deep Systemic Reset)..... | 28 |
| Usage Recommendations for Program Mode..... | 28 |

| | |
|--------------------------------------|----|
| 30 Factory Reset..... | 29 |
| Women’s Health..... | 29 |
| 33 UTI Kidney..... | 29 |
| 34 Skin Wound | 29 |
| 35 Anemia Care..... | 29 |
| 36 Cervix HPV | 30 |
| 37 Menstrual Care..... | 30 |
| 38 Vaginitis..... | 30 |
| 39 PID Adnexitis..... | 30 |
| 40 Urinary Incont..... | 30 |
| 41 STIs RTIs..... | 30 |
| 42 Ovarian Cysts..... | 30 |
| 43 Menopause Help..... | 30 |
| Other Specialized Health Modes | 31 |
| 44 Thyroid Care | 31 |
| 45 Bone Health | 31 |
| 46 GallBldr Pain | 31 |
| 47 Migraine Relief | 31 |
| 48 Sinus/Vertigo Focus | 31 |
| 49 Eyes Support..... | 31 |
| 50 Liver Support..... | 31 |
| Diseases and Frequency Set | 32 |
| Product Warranty (Hardware)..... | 46 |
| Limited Hardware Warranty | 46 |
| Warranty Claim Process: | 47 |
| Disclaimer of Medical Warranty | 47 |
| Limitation of Liability | 48 |
| Safety Disclaimer (Critical) | 48 |

Software/Firmware Disclaimer49
No Endorsement Disclaimer 50

Discover the Power of Zenthera Zapper



Zenthera Zapper brings ancient healing wisdom into modern technology with **50 precision-programmed frequency modes** that target everything from common infections and chronic pain to energy balancing and pet wellness. Inspired by Dr. Royal Rife's frequency therapy, Dr. Hulda Clark's protocols, and Dr. Piyush Saxena's cleansing sequences, this compact device delivers targeted electrical frequencies (Rife modes) or soothing sound vibrations (Sound Healing modes) through simple electrode or External Speaker with AUX Cable connections. Whether you're seeking daily detoxification, viral support, joint pain relief, or deep relaxation, Zenthera Zapper automates complete treatment protocols with crystal-clear LCD feedback showing real-time frequency, countdown timer, and progress.

Zenthera Zapper is safe for Humans / Pets and Food, it runs on 5v to 6v DC current. Which will not feel until you use it on thin skin like stomach.

No guesswork, no complex setup – just press MODE and heal. With intuitive 3-button operation, auto-advancing 7-minute sessions, and smart features like frequency skipping (UP/DOWN) and session memory, Zenthera Zapper fits seamlessly into your wellness routine. Run the **28-minute** Program Mode daily for systemic cleansing, choose from a wide range of electrode-based health protocols (roughly 21–240 minutes, depending on the mode), or relax with Sound Healing modes (about 5 -112 minutes) featuring Chakra tones, Brainwave entrainment, and Solfeggio-style frequencies. Safe, non-

invasive, and complementary to medical care, Zenthera Zapper empowers you to take charge of your natural healing journey – one frequency at a time.

Safety Guidelines

- **Use as directed** - complementary to medical care
- **3 / 7-minute sessions** prevent fatigue
- **Rest 20+ minutes** between modes
- **Hydrate well** during use
- **NOT for:** Pacemakers / pregnancy (consult doctor)
Stop if: Dizziness, discomfort occurs
- **Medical Disclaimer:** Zenthera Zapper provides experimental frequencies. This device provides experimental wellness frequencies; **no guaranteed results**; and not medical device.
- **Implants/medical devices:** Do not use if you have a pacemaker/ICD, implanted neurostimulator, insulin pump, cochlear implant, or any implanted electronic/metal device unless your doctor approves.
- **Heart/neck placement:** Never place electrodes across the chest/heart (left hand to right hand) or on the front of the neck.
- **Children & supervision:** Not for children without adult supervision/doctor guidance.
- **Skin safety:** Do not use on broken skin/open wounds with electrodes; stop if burning/tingling/pain occurs; clean electrodes after use.
- **Driving/sleep warning for sound modes:** Do not use Sleep/Brainwave sound healing modes while driving / operating machinery.
- **No replacement for treatment:** Strengthen to: Do not delay emergency care (chest pain, breathlessness, high fever, heavy bleeding, severe abdominal pain).
- **If symptoms worsen:** Stop using and consult a doctor if symptoms worsen or persist >48 hours.
- **Charging while using:** Do not use the device while charging.
- **Heat/fire:** Do not leave charging unattended, do not charge near pillow/blanket, stop charging if battery/device gets hot or swollen.

Battery Care Instructions

Power Source: 7.2V 1000 mAh internal Li-ion battery (about 7.2 Watt-hours).

Airline Travel Precautions

1. **Check airline policy** before flying. Most allow removable batteries <100Wh (Watt-hours).
2. **Remove battery for security:**
 - Unscrew 4 screws (keep loose, don't tighten)
 - Disconnect battery clip, remove battery
 - Place **Zapper body + loose screws** in one transparent zip bag
 - Place **Battery** in separate transparent zip bag
3. **No screwdriver** in carry-on/handbag
4. **Security demo:** Show how battery installs if asked

Battery Status Indicators

- LCD contrast low → Battery ~10%
- LCD barely visible → Battery ~5%
- LCD completely off → Battery 0% (dead)

Note: The zapper shows “**Batt Left: XX%**” at boot and has a critical low-battery lockout. when battery is repeatedly detected below the required power, the device beeps and displays “**Battery LOW / CHARGE NOW**”, stops zapping , and blocks until battery rises to minimum required power, then returns to menu.

if battery is too low, device will show “**Battery LOW! / CHARGE NOW**”, stop output, and won't run until battery recovers (the code blocks until voltage rises above a threshold, then returns to menu).

Charging Instructions

1. **Power OFF** device completely
2. Connect **9V 0.5A-1A AC-DC adapter** to power jack (use a good-quality 9V adapter (stable output). Avoid cheap unbranded adapters.)
3. Charge for **~2.0 hours maximum - avoid overcharging**
4. Disconnect charger when done


Pro Tips

- Charge when LCD contrast fades
- Keep spare battery for long sessions
- Store at 50-60% charge for longevity
- Charge every 2-3 months if unused

User Operating Manual

Welcome to Zenthera Zapper

Zenthera Zapper delivers precise healing frequencies through electrodes (health modes) or sound vibrations (Sound Healing modes). With **49 specialized modes**, it covers everything from infections and pain relief to chakra balancing and pet care. Each mode runs automatically with a clear LCD display showing frequency, time remaining, and progress.



Simple 3-Button Operation:

- **MODE button:** Select modes / Start-Stop sessions

- **UP button:** Previous mode / Next frequency
- **DOWN button:** Next mode / Previous frequency

Quick Start Guide

1. Power On

"Zenthera Zapper" and its firmware details ← Shows for 2 seconds

"DO NOT USE WITH PACEMKR/PREGNNCY" ← Warning will be displayed for 2.5 seconds.

"Batt Left: ~XX%" ← Battery screen for 1.5 seconds

"Select Mode: 01 Program Mode" ← Automatically shows.

Press **UP** or **DOWN** to change modes, **MODE** to start session

2. Select Mode

1. Press **MODE** → Shows menu: Select Mode: *"01 Program Mode"*
2. Use **UP/DOWN** to browse different modes
3. Last used mode remembers position

3. Start Session

Press **MODE** again → Session begins automatically:

Vir: 000.728 kHz ← Mode + Current Frequency

Tim:07:23 Rem:13 ← Time remaining + Frequencies left

4. During Session

- **UP:** Skip to next frequency (resets 7-min timer)
- **DOWN:** Go back to previous frequency
- **MODE:** Stop session (long beep)

Display Explained

Com: 000.727 kHz ← "Com" = Common Toxins mode, 727 Hz

Tim:06:45 Rem:12 ← 6 min 45 sec left, 12 frequencies remaining

Chk: 432 Hz OM ← Chakra mode showing "OM" label

Tim:07:12 Rem:0 ← Last frequency in set

- **Tim:** Remaining time for current frequency
- **Rem:** Frequencies left in complete set (most modes)



Using Carbon Electrodes (Pads / Straps)

- Use the supplied **carbon rubber electrodes** for all Health Modes that require electrodes (01-18, 02, 23-26, 31-47). The Sound Healing Modes 19–22, 27, 28 use an external speaker with AUX cable, not electrodes.
- Each pad has a **black carbon side** and a **blue plastic (insulated) side**. The black carbon side must touch the skin; the blue plastic side always faces outward.

- Apply a thin, even layer of conductive gel only on the **black carbon side** before placing it on the skin. Do not apply gel on the blue plastic side.
- Skin should be clean and dry (no oil, lotion, or broken skin). Do not place pads on cuts, rashes, infected or irritated areas. If you feel strong burning, sharp pain, or unusual discomfort under any pad, stop the session immediately and adjust pad position, gel, or end the session.
- For general whole-body support (Program, Common Toxins, Parasite Care, Detox Body, Immune Boost, etc.), place one pad on each wrist or each ankle so that the current passes through the body between the two pads.
- For local support (knee, shoulder, back, specific joint or area), place both pads so that the line between them “passes through” the problem area (for example, either side of the knee, above and below a painful spot on the spine, or front and back of a shoulder).
- Do not place electrodes **across the chest** (for example, one pad on the left chest and one on the right chest), on the front of the neck, or on the head/face.
- Use only on intact skin. Do not use if you have a pacemaker or are pregnant; consult your doctor before use for any serious medical Condition. Stop use if you experience dizziness, chest discomfort, breathing difficulty, or any unusual symptoms.
- After each session, gently remove the pads, wipe off excess gel from the skin and pad surface with a slightly damp cloth, squeeze out extra water from the pads if you rinsed them, and let the pads air-dry before storing. Ensure the device output jacks and cables stay dry at all times.

Core Modes Summary

Daily Essentials (Use Electrodes)

| Mode | Purpose | Total Time | Best For |
|-------------------|----------------------------|----------------|----------------|
| 01 Program | Dr. Piyush Saxena Protocol | 28 min | Daily cleanse |
| 04 Common | Cold, fever, infection | 111 min | General health |
| 05 Viral | Flu, COVID, dengue | 98 min | Viral support |
| 23 Detox | Body cleansing | 35 min | Weekly detox |

Specialized Health (Use Electrodes)

| Mode | Target | Frequencies | Total Time |
|------------------------|--------------------------------------|-------------|------------|
| 06 ENT | Ear, sinus, throat | 12 | 84 min |
| 07 Dental | Toothache, gums | 13 | 91 min |
| 08 Cardiac | Heart, circulation | 11 | 77 min |
| 13 Arthritis | Joint pain, back | 12 | 84 min |
| 14 Parasite | Worms, malaria | 23 | 161 min |
| 29 Rife Classic | Classic Frequencies of Rife Protocol | 11 | 77 min |
| 02 Core Support | Complete Core Frequency Set | 121 | 242 min |

Additional Health Modes (Use Electrodes)

| Mode | Target | Frequencies | Total Time |
|--------------------------|---|-------------|------------|
| 32 Pemphigus Vul | Autoimmune skin blistering support | 14 | 98 min |
| 33 UTI Kidney | UTI / E. coli urinary support | 36 | 108 min |
| 34 Skin Wound | Skin/wound infection support | 10 | 70 min |
| 35 Anemia Care | Anaemia / blood support set | 43 | 301 min |
| 36 Cervix HPV | Cervical HPV support | 16 | 112 min |
| 37 Menstrual Care | Menstrual/endometriosis/fibroid support | 18 | 126 min |
| 38 Vaginitis | Vaginitis/vaginosis/Candida support | 22 | 154 min |
| 39 PID Adnexitis | PID/adnexitis support | 7 | 49 min |

| Mode | Target | Frequencies | Total Time |
|---------------------------|---|--------------------|-------------------|
| 40 Urinary Incont | Urinary control/bladder support | 11 | 77 min |
| 41 STIs RTIs | STI/RTI support set | 15 | 105 min |
| 42 Ovarian Cysts | Ovarian cyst support | 3 | 21 min |
| 43 Menopause Help | Menopause support | 16 | 112 min |
| 44 Thyroid Care | Thyroid support | 11 | 77 min |
| 45 Bone Health | Bone regeneration support | 16 | 112 min |
| 46 GallBlDr Pain | Gallbladder support | 41 | 123 min |
| 47 Migraine Relief | Migraine / Headache relief | 5 | 35 min |
| 48 Sinus/Vertigo | Sinus, Sinusitis, Rhinitis, Tinnitus, Vertigo Support | 93 | 186 min |
| 49 Eyes Support | Cataracts, Glaucoma, Macular Degeneration, Eye Strain, Dry Eyes, Diabetic Retinopathy, Conjunctivitis, Blurred vision | 63 | 186 min |
| 50 Liver Support | Detox Of Anesthesia, Fatty Liver, Hangover, Jaundice, Lipoma, Hepatitis (A/B/C/Chronic/General), Listerioses, Liver Cirrhosis, Liver Diseases, Liver Enlargement, Parasites Flukes Liver, Pesticide Detox, Vitiligo, Yellow Fever | 116 | 348 min |

Sound Healing Modes (Use External Speaker with AUX Cable - No Electrodes)

| Mode | Purpose | Total Time |
|-----------------------|-------------------------|-------------------|
| 19 7-Chakras | Energy balancing | 42 min |
| 20 Beej Mantra | Subtle energy | 49 min |
| 21 Brain Wave | Sleep, focus | 98 min |
| 22 Healing | Pain relief, DNA repair | 112 min |

| | | |
|-----------------------------|------------------------------|----------|
| 27 Pets & Food | Pets/water energizing | 5-30 min |
| 28 Sleep & Relax | Helps to Sleep Stress relief | 55 min |

Note: Many frequencies (2-45 Hz) are too low to hear but work through vibration.

Special Modes

03 Sweep Mode (60 min)

- **Automatic scan:** 2 Hz → 500 kHz
- Covers full healing spectrum

26 ParasiteSweep (30 min)

- **Specialized scan:** 20 Hz → 50 kHz
- Targets parasite frequency range

31 Manual Mode (7 min)

- Start: 30 kHz
- **UP:** Increases by 1 Hz below 1000 Hz, otherwise by 10 Hz
- **DOWN:** Decreases the same way
- Perfect for custom frequencies, if you want to reset the frequency back to 30 kHz, then follow the **Factory Rest** process.

Button Functions by State

| State | MODE | UP | DOWN |
|---------|------------------|---------------|---------------|
| Idle | Enter Menu | - | - |
| Menu | Start Session | Previous Mode | Next Mode |
| Running | STOP (long beep) | Next Freq | Previous Freq |

Complete Session Flow Example

1. Power ON → "Press MODE to Start"
2. MODE → "Select Mode: 01 Program Mode"
3. UP/DOWN → Find desired mode
4. MODE → "Prg: 002.128 kHz" / "Tim:07:00 Rem:3"
5. Waits 7 min → Auto advances: "Prg: 10.000 kHz"
6. UP anytime → Skip to next frequency immediately
7. After 4th frequency → "Program Done!" → Back to Start

Daily Usage Recommendations

Beginner (30 min/day)

Day 1: 01 Program Mode (28 min)

Day 2: 04 Common Toxins (run only first 4 steps ≈ 28 min; use **UP** to skip to the next frequency)

Day 3: 23 Detox Body (35 min)

Note: If you feel discomfort, press MODE to stop and resume next day

Regular / Maintenance (45–120 min/day)

Option A (short): 01 Program (28 min) + one Health Mode for 15–60 min (use **UP** to shorten long sessions).

Option B (full): 01 Program (28 min) + one full Health Mode session (typical 60–120 min depending on mode).

Note: Long modes are optional to finish in one sitting; you can stop anytime with MODE and continue later.

Intensive Protocol

Morning: 01 Program (28 min) + 14 Parasite Care (shortened using **UP**, e.g., run 4–8 steps first; do not try to complete the whole mode in one go initially).

Evening: 04 Common Toxins (shortened if needed using **UP**) + any preferred Sound Healing mode (19–22, 28) using external speaker (no electrodes).

Note: Rest 20 minutes between long sessions and hydrate well.

External Speaker with AUX Cable and Volume Control Usage (Sound Healing Modes 19-22, 27, 28)

Disconnect electrodes completely. Connect External Speaker with AUX Cable/speaker. These gentle frequencies (2-963 Hz) work through calming sound vibrations. Many are subsonic (too low to hear) but powerful for relaxation, chakra balance, and pet calming. Use headphones for personal sessions or speakers for room therapy.

| Mode | Volume | Speaker Type | Best Setup |
|--------------------|-------------|--|--------------------------------|
| 19-20 Chakras/Beej | Low-Medium | Headphones | Personal energy work |
| 21 Brain Wave | Medium | Stereo speakers | Binaural entrainment |
| 22 Healing | Medium-Loud | Room speakers | Solfeggio immersion |
| 27 Pets & Food | Loud | Bowl next to speaker | Vibrational transfer |
| 28 Sleep & Relax | Low-Medium | Headphones or Stereo speakers near the bed | Restful slumber, stress relief |

Selecting the external Speakers for Sound Healing Modes

External speakers for Zenthera Zapper's Sound Healing modes (19–22, 27& 28) should be chosen to reproduce 2 Hz – 1 kHz reliably at comfortable volume, using the External Speaker with AUX Cable (no electrodes).



Must-have specifications

- **Input:** 3.5 mm AUX/line-in (stereo preferred).
- **Speaker type:** Active (powered) speakers or any speaker connected to a small amplifier (phone-passive speakers alone won't work).
- **Frequency response** (realistic): At least 40 Hz - 10 kHz on paper; lower is better, but most consumer speakers won't truly reproduce 2–20 Hz (those will be felt as vibration more than heard).
- **Low noise floor:** minimal hiss at low volume, since these modes are used at low–medium volume and often during relaxation/sleep.
- **Volume Control:** The Speaker must have volume control.

Recommended (best results) specs

- **Stereo playback:** For Brain Wave (Mode 21) and Sleep/Relax-style entrainment, use stereo speakers (L/R) or stereo headphones, because the intended setup for brainwave entrainment is stereo.
- **Bass capability:** a speaker system that can reach ~40 Hz (or a 2.0 + subwoofer system reaching ~30–35 Hz) helps you feel the low-frequency components more clearly.
- **Power:** roughly 10–50 W RMS total is enough for a bedroom/room at low to medium volume without distortion.

What to avoid

- Bluetooth-only speakers with no AUX input, because the document's recommended connection method is the AUX cable.
- Very small/mono speakers if you want the full "Brain Wave" style effect; stereo is preferred for those modes.

Troubleshooting

- No display? → Battery Drained, turn the on/off button to OFF and connect the 9v AC-DC Adapter. Wait for 90 min.
- No sound/beeps? → Device working normally
- Session won't start? → Press MODE twice
- Frequency stuck? → Use UP or DOWN button

Support: Last mode auto-saves. Power cycle returns to your favourite mode.

Frequency Table by Mode

Note: Most modes run 7 minutes per step; some modes run **3 minutes** per step (e.g., Common Toxins, Parasite Care, Cancer Support, Rife Classics, UTI+Kidney, Anemia, Gallbladder, Eyes Support), and some run **2 minutes** per step (Core Support, Sinus/Vertigo).

| Mode | Illnesses Targeted | Rationale / Explanation | Total Time in minutes | Auto - Advance | Short Form |
|-------------------------|--|---|-----------------------|----------------|------------|
| 01 Program Mode | General Protocol (Dr. Piyush Saxena) | High-frequency protocol steps; 2128 common in systemic cleansing, 30kHz from Dr. Clark/Hulda Clark protocols | 28 | Yes | Prg |
| 02 Core Support | Broad-spectrum helper set | Complete backbone frequency set. 2 Minutes per frequency, total 121 frequencies. | 242 | Yes | CoS |
| 03 Sweep Mode | Broad-Spectrum Pathogen Scan | Full spectrum sweep technique used in Rife experimentation to hit unknown pathogens across entire frequency range | 60 | Continuous | Swp |
| 04 Common Toxins | Cough, Cold, Fever, Headache, Sore Throat, Allergy, Asthma, Bronchitis, Sinusitis, Diarrhea, Nausea, Fungus, Boils/Pimples (Skin Infection), General Infection | 37 frequencies running, 3 minutes each frequency. | 111 | Yes | Com |
| 05 Viral Relief | COVID, Dengue, Flu, Malaria, Hepatitis, Typhoid, Gastroenteritis | Classic antiviral trio + high harmonics for lung/viral load | 98 | Yes | Vir |
| 06 ENT Care | Ear Infection, Tinnitus, Sinus, Sore Throat, Tonsillitis, Vertigo | ENT infection core and tinnitus/hearing support | 84 | Yes | ENT |
| 07 Dental Care | Toothache, Gums, Abscess, Oral Thrush, Jaw Pain | Dental abscess protocol (465-880 core); 10000 for chronic infection | 91 | Yes | Den |

| Mode | Illnesses Targeted | Rationale / Explanation | Total Time in minutes | Auto - Advance | Short Form |
|----------------------------------|--|--|------------------------------|-----------------------|-------------------|
| 08 Cardiac Care | Arrhythmia, Angina, Blood Pressure, Circulation | Arrhythmia-specific (225/510); vascular support harmonics | 77 | Yes | Car |
| 09 Metabolic Care | Diabetes, Obesity, Thyroid, Kidney, Liver | Infection + liver/kidney detox; terrain support for metabolic stress | 98 | Yes | Met |
| 10 Organ Health (General) | Lungs, Stomach, Brain, Eyes, Skin, Pancreas, COPD | Multi-organ infection + detox; lung/brain inflammation. 3 min each frequency | 126 | Yes | Org |
| 11 Reproductive Organs | PCOS, Fibroids, UTI, HPV, Prostate, Infertility | Pelvic infection core; structural support for fibroids/cysts | 112 | Yes | Rep |
| 12 Autoimmune | Lupus, Psoriasis, Crohn's, Vitiligo, Celiac | Immune modulation + infection terrain; not disease-specific | 84 | Yes | Alm |
| 13 Arthritis Care | Joint/Back Pain, Sciatica, Fibromyalgia, Gout, OA/RA | Pain relief + arthritis core | 84 | Yes | Art |
| 14 Parasite Care | Worms, Giardia, Amoeba, Tapeworm, Malaria Parasite | Low fundamentals + high parasite harmonics (tapeworm range) – 3 min each frequency | 99 | Yes | Par |
| 15 Cancer Care | Breast, Lung, Prostate, Colon, Lymphoma Support | Cancer terrain + infection/pain support - 3 min each frequency | 84 | Yes | Can |
| 16 Ulcer Care | Peptic Ulcers, Stomach Ulcers, Fistula support (Dental fistula, Tracheoesophageal fistula, Fistula general), General Ulcer | H. pylori/infection core + stomach lining support; includes fistula/dental fistula support frequencies | 81 | Yes | Ulc |
| 17 Piles & Gut | Haemorrhoids, Gut Inflammation, Colitis | Colon/rectum infection + circulation support | 63 | Yes | P&G |
| 18 Paralysis | Nerve Recovery, Paralysis Support | Brainwave range + neuropathy inflammation support | 56 | Yes | Prl |
| 19 7-Chakras | Energetic balance, stress, emotional state (non-medical sound healing) | Chakra tone mappings from sound-therapy traditions; intended for relaxation and | 42 | Yes | Chk |

| Mode | Illnesses Targeted | Rationale / Explanation | Total Time in minutes | Auto - Advance | Short Form |
|---------------------------|---|---|------------------------------|-----------------------|-------------------|
| | | balance, not a medical frequency set. | | | |
| 20 Beej Mantra | Subtle energy, mantric resonance (non-medical sound healing). | Seed-syllable (“Beej”) tone mappings; described as experimental/spiritual heuristic practice / treatment. | 49 | Yes | BJM |
| 21 Brain Wave | Sleep, anxiety, cognition, emotional balance (sound healing). | Uses brainwave ranges (Delta/Theta/Alpha/Beta/Gamma) and Solfeggio-style tones commonly used in binaural/entrainment discussions for calming and focus. | 98 | Yes | Brn |
| 22 Healing Mode | Pain relief, tissue healing, relaxation, emotional healing (sound healing). | Solfeggio-tradition “healing tones” bundle for mood/relaxation and supportive pain/tissue themes; positioned as non-medical sound healing. | 112 | Yes | Hea |
| 23 Detox Body | General Detoxification, Cleansing | Solfeggio cleanse + Rife detox core | 35 | Yes | Det |
| 24 Immune Boost | Immunity, Lymphatic Drainage, Circulation | Lymphatic/immune helper frequencies | 35 | Yes | Imm |
| 25 Hormonal Care | Thyroid, Adrenal, Reproductive Hormones | Endocrine infection support + speculative hormone tones | 70 | Yes | Hor |
| 26 Parasite Sweep | Broad Parasite Spectrum (Worms, Flukes, Protozoa) | Targeted parasite range sweep; covers common parasite harmonics | 30 | Continuous | PSw |
| 27 Pets & Food | Pet Support, Food/Water Energizing (sound healing) | Broad-spectrum frequency for pet wellness and food/water vibrational treatment UP/DOWN selects duration (default 15 min) | 5, 10, 15, 20, 25, 30 min | Fixed | P&F |

| Mode | Illnesses Targeted | Rationale / Explanation | Total Time in minutes | Auto - Advance | Short Form |
|-----------------------------|---|---|------------------------------|-----------------------|-------------------|
| 28 Sleep & Relax | Stress relief, Restful slumber (sound healing) | Progressive brainwave entrainment: Alpha → Delta + deep relaxation | 55 | Yes | S&R |
| 29 Rife Classics | BX virus, staph/strep, foundational pathogens | Original Rife lab frequencies - 3 min each frequency | 70 | Yes | Rif |
| 30 Factory Reset | Erases flash memory (menu + frequency) | Clears saved settings (last mode + saved manual frequency) and restarts the device. | N/A | N/A | FRE |
| 31 Manual Mode | Custom Single Frequency | User-defined; with default 30kHz is classic Hulda Clark zapper frequency for broad-spectrum pathogen disruption | 7 | Fixed | Man |
| 32 Pemphigus Vul | Autoimmune skin blistering support (Pemphigus vulgaris support) | Condition-focused supportive set intended for autoimmune/skin healing support (complementary). | 98 | Yes | PeV |
| 33 UTI Kidney | UTI / E. coli urinary support, kidney/urinary discomfort support. | Focused urinary-tract infection support set, meant to simplify selection vs general modes. | 108 | Yes | UTI |
| 34 Skin Wound | Skin/wound infection support. | Condition-focused infection + skin healing support set aimed at common wound/skin infection patterns. | 70 | Yes | SKN |
| 35 Anemia Care | Anemia / blood support. | Long supportive set intended for blood/iron/anemia-type support themes (complementary). | 129 | Yes | ANM |
| 36 Cervix HPV | Cervical HPV support. | Focused cervix/HPV supportive set designed for quicker selection than general reproductive mode. | 112 | Yes | HPV |

| Mode | Illnesses Targeted | Rationale / Explanation | Total Time in minutes | Auto - Advance | Short Form |
|--------------------------|--|---|------------------------------|-----------------------|-------------------|
| 37 Menstrual Care | Menstrual issues, endometriosis/fibroid support. | Reproductive / uterine supportive frequency group tailored to menstrual irregularity and fibroid / endometriosis-type support themes. | 126 | Yes | MNC |
| 38 Vaginitis | Vaginitis/vaginosis/Candida support. | Focused vaginal flora / yeast - support frequency grouping; positioned as supportive. | 154 | Yes | VGN |
| 39 PID Adnexitis | PID/adnexitis support. | Short targeted pelvic inflammatory support set for fast selection. | 49 | Yes | PID |
| 40 Urinary Incont | Urinary control / bladder support. | Supportive bladder/urinary control frequency grouping (symptom / support oriented). | 77 | Yes | UIC |
| 41 STIs RTIs | STI/RTI support set (sexual / reproductive tract infection support). | Focused reproductive infection support grouping intended to cover common STI / RTI-type support patterns (complementary). | 105 | Yes | STI |
| 42 Ovarian Cysts | Ovarian cyst support. | Very short targeted ovarian-cyst supportive set for quick use. | 21 | Yes | OVC |
| 43 Menopause Help | Menopause support. | Symptom/support-oriented menopause frequency grouping (complementary). | 112 | Yes | MNP |
| 44 Thyroid Care | Thyroid support. | Condition-focused thyroid supportive set (not a replacement for medical thyroid treatment). | 77 | Yes | THY |
| 45 Bone Health | Bone regeneration / bone health support. | Supportive bone healing/regeneration frequency set for bone health themes. | 112 | Yes | BON |

| Mode | Illnesses Targeted | Rationale / Explanation | Total Time in minutes | Auto - Advance | Short Form |
|---------------------------|---|---|------------------------------|-----------------------|-------------------|
| 46 GallBldr Pain | Gallbladder Pain Support | Condition-focused gallbladder supportive grouping intended for gallbladder discomfort themes. | 123 | Yes | GAL |
| 47 Migraine Relief | Migraine / Headache relief | Headache support | 35 | Yes | MIG |
| 48 Sinus / Vertigo | Sinusitis (acute/chronic), rhinitis/runny nose, tinnitus (ringing), vertigo/giddiness, ear pressure/ENT congestion themes. | Focused ENT subset combining sinus + nasal + inner-ear support frequencies commonly repeated across broader sets, intended as a more targeted option than general ENT when symptoms are mainly sinus blockage, nasal irritation, ringing, or dizziness. | 186 | Yes | SRV |
| 49 Eyes Support | Cataracts, Glaucoma, Macular Degeneration, Eye Strain, Dry Eyes, Diabetic Retinopathy, Conjunctivitis, Blurred vision | Eyes Support is an frequency set intended to provide general eye-comfort support (surface irritation/strain) and broad eye-health support themes (lens/retina/optic nerve) | 186 | Yes | Eye |
| 50 Liver Support | Hangover, Jaundice, Lipoma, Listerioses, Liver stores/bile/cholesterol support, Liver enlargement, Hepatitis A/B/C (including chronic, cryptogenic, viral), Liver cirrhosis, Liver diseases, Fatty liver, Parasites flukes liver, Pesticide detox, Surgery anesthesia detox, Vitiligo, Yellow fever | Comprehensive liver health and detoxification support covering multiple hepatitis types, cirrhosis, fatty liver, jaundice, hangover recovery, pesticide detoxification, parasitic liver flukes, and post-anesthesia recovery. Complementary to medical treatment, not a replacement. For acute liver conditions, seek immediate medical care. | 116 | Yes | Liv |

Detailed Sweep & Manual Mode Explanations

Mode 03: Sweep Mode

- **Formula:** 2 Hz → 500 kHz over 60 minutes
- "Full spectrum sweeps" are a **classic technique** for unknown pathogens.
- **Rationale:** Covers **all known pathogen ranges**:

| Range | Targets |
|-------------|--|
| 2-1000 Hz | Pain, Inflammation, Low-Frequency Bacteria |
| 1000-10kHz | Most Viruses, Common Bacteria, Candida |
| 10k-100kHz | Parasites, Flukes, Advanced Infections |
| 100k-500kHz | Experimental High-Range Pathogens |

Mode 26: Parasite Sweep

- **Formula:** ~20 Hz → 50 kHz over 30 minutes
- Parasite frequencies cluster between **400 Hz-50 kHz** (tapeworms 15-40kHz, flukes 20-30kHz).
- **Rationale:** More targeted than full sweep; focuses on documented parasite harmonic ranges while avoiding ultra-low/high extremes.

Mode 31: Manual Mode

- **Formula:** Starts at 30 kHz (Dr. Clark standard); changes by +1 Hz / -1 Hz below 1000 Hz, else +10 Hz / -10 Hz (2 Hz minimum, 500 kHz maximum)
- 30 kHz is the **Hulda Clark zapper frequency** – specifically chosen because many pathogens resonate/die at this square-wave frequency.
- **Rationale:** Allows experimentation with exact frequencies from Rife lists while maintaining 7-minute session safety timer. Perfect for targeting specific pathogens identified in research.

Special Modes (Non-Medical / Sound Healing)

Sound Healing Modes (19-22, 27, 28): For Chakras, Brain Waves, Healing, and Pets & Food, use the External Speaker (with volume control) with AUX Cable instead of electrodes.



These gentle frequencies work through sound vibrations. Many are too low to hear but still help relax your mind and body. Connect the External Speaker with AUX Cable, play through External Speaker with AUX Cable, and let the calming tones do their work.

| Mode | Focus | Source | Explanation | Total Duration in minutes | Auto-Advance |
|---------------------------|--|----------------------|---|---------------------------|--------------|
| 19 7-Chakras | Energy Center Balance | Sound healing | Chakra tone mappings | 42 | Yes |
| 20 Beej Mantra | Subtle Energy / Mantras | Spiritual heuristics | Seed syllable frequencies; experimental | 49 | Yes |
| 21 Brain Wave | Sleep, Anxiety, Cognition | Binaural + Solfeggio | Brainwave entrainment + emotional healing tones | 98 | Yes |
| 22 Healing Mode | Pain, Tissue Repair, Emotional Healing | Solfeggio tradition | Extended healing frequency set; mood/dopamine support | 112 | Yes |
| 27 Pets & Food | Pet Support, Food Energizing | General support | Broad-spectrum frequency for pet wellness and | 5 to 30 | Fixed |

| Mode | Focus | Source | Explanation | Total Duration in minutes | Auto-Advance |
|-----------------------------|--------------------------------|---------------|--|---------------------------|--------------|
| | | | food/water vibrational treatment | | |
| 28 Sleep & Relax | Stress relief, Restful slumber | Sound healing | Progressive brainwave entrainment: Alpha → Delta + deep relaxation | 55 min | Yes |



01 Program Mode

- **Dr. Piyush Saxena's Protocol:** These exact frequencies (2128 Hz, 10 kHz, 30 kHz, 48 kHz) form his core "Cleansing Therapy" sequence for general detoxification, pathogen clearance, and systemic reset. Dr. Saxena promotes this as a foundational 4-step program for natural healing without drugs, often run as 7-minute steps with breaks.
- **Standards:**
 - **2128 Hz:** Ubiquitous "mortal oscillatory rate" (MOR) for broad-spectrum pathogens (bacteria, viruses, abdominal infections).
 - **10 kHz:** General antimicrobial, detoxification support (Dr. Clark/Hulda Clark influence).
 - **30 kHz:** Dr. Hulda Clark's universal parasite/virus killer (7-20-7-20-7 protocol at 30 kHz).

- **48 kHz:** High-frequency systemic cleanse, cancer support, and deep pathogen targeting.

Protocol Structure: Each runs for 7 minutes (total ~28 minutes), auto-advancing. UP/DOWN skips steps. This matches empirical Rife usage (e.g., 45+ min on 2128 Hz for cancer) but is optimized for daily use.

Scientific Context: Rife therapy targets microbial resonance via electromagnetic fields. While anecdotal/clinical, no mainstream medical validation exists. Cancer Research UK notes no cancer cure evidence. Use as complementary only.

Detailed Illness Associations and Frequency Explanations

Step 1: 2128 Hz (Dr. Piyush Saxena's Primary Frequency)

Targeted Illnesses/Conditions:

- **General Infections:** Abdominal infections, chlamydia, cytomegalovirus (CMV), ear disorders, intestinal disorders.
- **Cancer Support:** Breast cancer (primary/secondary), prostate cancer, sarcoma, melanoma (use 56+ min).
- **Prostatitis/BPH:** 10+ min essential.
- **Chronic Conditions:** Allergies, E. coli, hookworm, hypertension support.

Mechanism: Rife's "tertiary cancer frequency" and broad-spectrum killer. Disrupts pathogen cell walls via resonance. Dr. Saxena starts here for foundational cleansing.

Step 2: 10 kHz (Systemic Antimicrobial)

Targeted Illnesses/Conditions:

- **Athlete's Foot, Aura Cleansing:** Fungal/mold support.
- **Dental Infections, Dandruff:** Secondary oral/skin pathogens.
- **General Detox:** Bacteria, prostatitis, colorectal issues.
- **Cancer Adjunct:** Prostate, general carcinoma.

Mechanism: High-energy frequency for biofilm disruption and deep-tissue penetration. Complements 2128 Hz for multi-layered pathogen elimination.

Step 3: 30 kHz (Dr. Hulda Clark Universal Killer)

Targeted Illnesses/Conditions:

- **Parasites/Viruses:** All parasites (Worms, Flukes, Protozoa), viruses (EBV, HIV support), Bacteria Generally.
- **Liver/Gallbladder Cleanse:** Dr. Saxena's heart attack prevention protocol post-cleanse.
- **Immune Boost:** Typhoid, Klebsiella, General infections.

Mechanism: Hulda Clark's single-frequency protocol (30 kHz for 7 min x3 with 20-min breaks) destroys all Pathogens via Square-Wave Zapping. Simplifies sets; Dr. Saxena adapts for daily use.

Step 4: 48 kHz (Deep Systemic Reset)

Targeted Illnesses/Conditions:

- **Cancer (Advanced):** BX cancer virus (use with 21275 Hz if capable), lymphoma support.
- **High-Frequency Pathogens:** Resistant strains, deep-tissue infections not hit by lower freqs.
- **Prostate/Colorectal:** Final cleanse phase.

Mechanism: Ultra-high frequency for "finishing" resistant microbes. Targets high-resonance pathogens missed by lower bands. Completes Saxena's protocol for full-spectrum coverage.

Usage Recommendations for Program Mode

- **Daily:** Run full 28-min cycle 1-2x/day for maintenance.
- **Intensive:** Extend 2128 Hz (UP to repeat) for cancer/chronic issues.
- **Integration:** Follow with Common Toxins or Parasite Care modes.
- **Disclaimer:** Experimental/complementary. Not FDA-approved; consult physician. No proven cure for diseases like cancer.

This validation confirms the frequencies align with Dr. Saxena and Rife-style frequency traditions, based on multiple independent references and commonly used frequency lists.

30 Factory Reset

Focus: If you want to reset the Zenthera Zapper to restore Manual Mode back to 30 kHz and clear the last-used menu selection, use the Factory Reset option. This erases saved last mode + last manual frequency and restarts the device.

Steps:

- From the mode selection menu, navigate to “**30 Factory Reset**” and press MODE to enter the reset confirmation screen.
- The device will display: “**Factory Reset?**” and “**UP:Y DOWN/MODE:N**”.
- Press UP button to confirm (triple-beep confirmation) and perform the Reset + Reboot, or press DOWN or MODE button to Cancel and return to the menu.
- Factory reset does not change hardware, only clears memory/settings.

Women’s Health

Note: In addition to broad modes, the firmware includes dedicated women’s health and Condition-focused modes for faster selection:

33 UTI Kidney

Illnesses/Focus: UTI / urinary tract support (E. coli–focused set).

34 Skin Wound

Illnesses/Focus: Skin/wound infection support.

35 Anemia Care

Illnesses/Focus: Anemia / blood support.

36 Cervix HPV

Illnesses/Focus: Cervical HPV support.

37 Menstrual Care

Illnesses/Focus: Menstrual disorders; endometriosis and fibroid/cyst support.

38 Vaginitis

Illnesses/Focus: Vaginitis/vaginosis; Candida/yeast support.

39 PID Adnexitis

Illnesses/Focus: PID (pelvic inflammatory disease) / adnexitis support.

40 Urinary Incont

Illnesses/Focus: Urinary incontinence; bladder/prostate/urinary control support.

41 STIs RTIs

Illnesses/Focus: STI/RTI support (includes gonorrhea/chlamydia/cervicitis-oriented set).

42 Ovarian Cysts

Illnesses/Focus: Ovarian cyst support.

43 Menopause Help

Illnesses/Focus: Menopause support.

Other Specialized Health Modes

44 Thyroid Care

Illnesses/Focus: Thyroid support (hypothyroidism, hyperthyroidism, thyroiditis).

45 Bone Health

Illnesses/Focus: Bone regeneration; fractures/osteoporosis/bone healing support.

46 GallBldr Pain

Illnesses/Focus: Gallbladder support (gallbladder pain/gallstones-oriented set).

47 Migraine Relief

Illnesses/Focus: Migraine / Headache relief.

48 Sinus/Vertigo Focus

Illnesses/Focus: Sinusitis (acute/chronic), rhinitis/runny nose, tinnitus (ringing), vertigo/giddiness, ear pressure/ENT congestion themes.

49 Eyes Support

Illnesses/Focus: Cataracts, Macular Degeneration, Glaucoma, Eye Strain, Dry Eyes, Diabetic Retinopathy, Conjunctivitis, Blurred vision

50 Liver Support

Illnesses/Focus: Comprehensive liver health and detoxification support covering multiple hepatitis types, cirrhosis, fatty liver, jaundice, hangover recovery, pesticide detoxification, parasitic liver flukes, and post-anesthesia recovery. Must run **14 Parasite Care along with 50 Liver Support.**

Diseases and Frequency Set

This table is a wellness-support lookup for choosing modes; it is not diagnostic, and results are not guaranteed. For emergencies or serious symptoms, seek medical care.

Note: Mode 02 Core Support has list of frequencies to support most of the Virus, Parasite, Cancer, Illness, Bacteria, Pathogen, Condition, Fungus etc.

| S/No | Disease/Pathogen | Type | Recommended mode(s) |
|------|-------------------------------------|-----------|---|
| 1 | Acanthamoeba infection | Parasite | 14 Parasite Care, 02 Core Support |
| 2 | Acinetobacter baumannii infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 3 | Acne | Bacteria | 34 Skin Wound, 04 Common Toxins |
| 4 | Actinomyces infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 5 | Acute bronchitis | Condition | 10 Organ Health, 04 Common Toxins/General Infection |
| 6 | Acute gastroenteritis | Condition | 16 Ulcer Care, 02 Core Support |
| 7 | Acute hepatitis | Condition | 02 Core Support |
| 8 | Acute lymphoblastic leukemia | Cancer | 15 Cancer Care, 02 Core Support |
| 9 | Acute myeloid leukemia | Cancer | 15 Cancer Care, 02 Core Support |
| 10 | Acute pancreatitis | Condition | 02 Core Support |
| 11 | Acute sinusitis | Condition | 48 Sinus/Vertigo, 06 ENT Care, 04 Common Toxins (if fever/infection), 02 Core Support |
| 12 | Acute viral hepatitis (unspecified) | Virus | 05 Viral Relief, 02 Core Support |
| 13 | acute viral respiratory infection | Virus | 05 Viral Relief, 02 Core Support |
| 14 | Adenovirus conjunctivitis | Virus | 05 Viral Relief, 02 Core Support |
| 15 | Adenovirus gastroenteritis | Virus | 05 Viral Relief, 02 Core Support, 49 Eyes Support |
| 16 | Adenovirus respiratory | Virus | 05 Viral Relief, 02 Core Support |
| 17 | Adrenal cancer | Cancer | 15 Cancer Care, 02 Core Support |
| 18 | Aeromonas infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 19 | Allergic rhinitis | Condition | 48 Sinus/Vertigo, 06 ENT Care, 02 Core Support |
| 20 | Amoeba liver abscess | Parasite | 14 Parasite Care, 02 Core Support |
| 21 | Amoebiasis | Parasite | 14 Parasite Care, 02 Core Support |

| S/No | Disease/Pathogen | Type | Recommended mode(s) |
|-------------|---|-------------|--|
| 22 | Anemia (iron deficiency) | Condition | 35 Anemia Care |
| 23 | Angina pectoris | Condition | 08 Cardiac Care, 02 Core Support |
| 24 | Anisakiasis | Parasite | 14 Parasite Care, 02 Core Support |
| 25 | Anthrax | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 26 | Anxiety disorder | Condition | 22 Healing Mode, 02 Core Support |
| 27 | Aphthous ulcers | Condition | 16 Ulcer Care, 02 Core Support |
| 28 | Arrhythmia | Condition | 08 Cardiac Care, 02 Core Support |
| 29 | Arthritis (osteo) | Condition | 13 Arthritis/Pain, 02 Core Support |
| 30 | Ascariasis | Parasite | 14 Parasite Care, 02 Core Support |
| 31 | Aspergillosis | Fungus | 04 Common Toxins/General Infection, 07 Dental Care |
| 32 | Asthma (chronic) | Condition | 04 Common Toxins/General Infection, 02 Core Support |
| 33 | Atrial fibrillation | Condition | 08 Cardiac Care, 02 Core Support |
| 34 | Avian influenza (H5N1 risk) | Virus | 05 Viral Relief, 02 Core Support |
| 35 | Babesiosis (rare) | Parasite | 14 Parasite Care, 02 Core Support |
| 36 | Bacillus cereus food poisoning | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 37 | Bacterial conjunctivitis | Bacteria | 04 Common Toxins/General Infection, 02 Core Support, 49 Eyes Support |
| 38 | Bacterial infection (bloodstream, unspecified) | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 39 | Bacterial infection (gastrointestinal, unspecified) | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 40 | Bacterial infection (post-surgical, unspecified) | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 41 | Bacterial infection (respiratory, unspecified) | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 42 | Bacterial infection (skin, unspecified) | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 43 | Bacterial infection (urinary, unspecified) | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 44 | Bacterial infection (wound, unspecified) | Bacteria | 34 Skin Wound, 02 Core Support |
| 45 | Bacterial meningitis | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 46 | Bacterial vaginosis | Bacteria | 38 Vaginitis, 02 Core Support |
| 47 | Bartonella henselae (cat-scratch) | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 48 | Basal cell carcinoma | Cancer | 15 Cancer Care, 02 Core Support |

| S/No | Disease/Pathogen | Type | Recommended mode(s) |
|------|-------------------------------------|-----------|---|
| 49 | Bedbug infestation | Parasite | 14 Parasite Care, 02 Core Support |
| 50 | Benign prostatic hyperplasia | Condition | 11 Reproductive Health, 02 Core Support |
| 51 | BK virus | Virus | 05 Viral Relief, 02 Core Support |
| 52 | Bladder cancer | Cancer | 15 Cancer Care, 02 Core Support |
| 53 | Blastomycosis (rare) | Fungus | 04 Common Toxins/General Infection, 02 Core Support |
| 54 | Blurred vision | Condition | 49 Eyes Support, 10 Organ Health, 02 Core Support |
| 55 | Boils/furunculosis | Bacteria | 34 Skin Wound |
| 56 | Bronchiectasis | Condition | 10 Organ Health, 04 Common Toxins/General Infection |
| 57 | Brucella melitensis infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 58 | Brucellosis | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 59 | Campylobacter enteritis | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 60 | Campylobacter jejuni infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 61 | Candidiasis | Fungus | 04 Common Toxins/General Infection, 07 Dental Care |
| 62 | Cataracts | Condition | 49 Eyes Support, 10 Organ Health, 02 Core Support |
| 63 | Cellulitis | Bacteria | 34 Skin Wound |
| 64 | Cervicitis / STI/RTI (general) | Condition | 41 STIs RTIs, 02 Core Support |
| 65 | Chandipura virus | Virus | 05 Viral Relief, 02 Core Support |
| 66 | Chikungunya | Virus | 05 Viral Relief, 02 Core Support |
| 67 | Chlamydia | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 68 | Chlamydomphila pneumoniae infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 69 | Cholangiocarcinoma | Cancer | 15 Cancer Care, 02 Core Support |
| 70 | Cholera | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 71 | Chronic diarrhea | Condition | 17 Piles/Gut, 02 Core Support |
| 72 | Chronic fatigue syndrome | Condition | 13 Arthritis/Pain, 22 Healing Mode |
| 73 | Chronic kidney disease | Condition | 09 Metabolic Care, 02 Core Support |
| 74 | Chronic lymphocytic leukemia | Cancer | 15 Cancer Care, 02 Core Support |
| 75 | Chronic myeloid leukemia | Cancer | 15 Cancer Care, 02 Core Support |
| 76 | Chronic pancreatitis | Condition | 02 Core Support |
| 77 | Chronic pelvic pain | Condition | 11 Reproductive Health, 25 Hormonal Care |
| 78 | Chronic prostatitis | Condition | 11 Reproductive Health, 25 Hormonal Care |
| 79 | Chronic sinusitis | Condition | 48 Sinus/Vertigo, 06 ENT Care, 02 Core Support |
| 80 | Chronic tonsillitis | Condition | 02 Core Support |

| S/No | Disease/Pathogen | Type | Recommended mode(s) |
|-------------|--|-------------|---|
| 81 | Cirrhosis | Condition | 50 Liver Support, 02 Core Support |
| 82 | Citrobacter freundii infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 83 | Clonorchiasis (rare) | Parasite | 14 Parasite Care, 02 Core Support |
| 84 | Clostridioides difficile infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 85 | Clostridium perfringens food poisoning | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 86 | CMV infection | Virus | 05 Viral Relief, 02 Core Support |
| 87 | Coccidioidomycosis (travel risk) | Fungus | 04 Common Toxins/General Infection, 02 Core Support |
| 88 | Community-acquired pneumonia | Condition | 10 Organ Health, 04 Common Toxins/General Infection |
| 89 | Constipation (functional) | Condition | 17 Piles/Gut, 02 Core Support |
| 90 | COPD | Condition | 10 Organ Health, 02 Core Support |
| 91 | Coronary artery disease | Condition | 08 Cardiac Care, 02 Core Support |
| 92 | Coronavirus 229E | Virus | 05 Viral Relief, 02 Core Support |
| 93 | Coronavirus HKU1 | Virus | 05 Viral Relief, 02 Core Support |
| 94 | Coronavirus NL63 | Virus | 05 Viral Relief, 02 Core Support |
| 95 | Coronavirus OC43 | Virus | 05 Viral Relief, 02 Core Support |
| 96 | COVID-19 | Virus | 05 Viral Relief, 02 Core Support |
| 97 | Coxiella burnetii (Q fever) | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 98 | Coxsackievirus A16 | Virus | 05 Viral Relief, 02 Core Support |
| 99 | Coxsackievirus B | Virus | 05 Viral Relief, 02 Core Support |
| 100 | Crimean-Congo hemorrhagic fever | Virus | 05 Viral Relief, 02 Core Support |
| 101 | Cryptococcosis | Fungus | 04 Common Toxins/General Infection, 07 Dental Care |
| 102 | Cryptosporidiosis | Parasite | 14 Parasite Care, 02 Core Support |
| 103 | Cutaneous larva migrans | Parasite | 14 Parasite Care, 02 Core Support |
| 104 | Cutaneous leishmaniasis | Parasite | 14 Parasite Care, 02 Core Support |
| 105 | Cyclospora infection | Parasite | 14 Parasite Care, 02 Core Support |
| 106 | Cystoisospora infection | Parasite | 14 Parasite Care, 02 Core Support |
| 107 | Dengue | Virus | 05 Viral Relief, 02 Core Support |
| 108 | Dental abscess | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 109 | Dental caries | Condition | 02 Core Support |
| 110 | Dental fistula (Fistula dentalis) | Condition | 16 Ulcer Care, 07 Dental Care (optional), 02 Core Support |
| 111 | Depression | Condition | 22 Healing Mode, 02 Core Support |
| 112 | Dermatitis/eczema | Condition | 10 Organ Health, 02 Core Support |

| S/No | Disease/Pathogen | Type | Recommended mode(s) |
|-------------|---|-------------|--|
| 113 | Dermatophytosis (ringworm) | Fungus | 04 Common Toxins/General Infection, 07 Dental Care |
| 114 | Diabetic neuropathy | Condition | 13 Arthritis/Pain, 22 Healing Mode |
| 115 | Diabetic Retinopathy | Condition | 49 Eyes Support, 09 Metabolic Care, 10 Organ Health, 02 Core Support |
| 116 | Diphtheria | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 117 | Dirofilariasis (rare) | Parasite | 14 Parasite Care, 02 Core Support |
| 118 | Dracunculiasis (eradicated but history) | Parasite | 14 Parasite Care, 02 Core Support |
| 119 | Dry eyes | Condition | 49 Eyes Support, 10 Organ Health, 02 Core Support |
| 120 | Dysmenorrhea | Condition | 11 Reproductive Health, 25 Hormonal Care |
| 121 | E. coli diarrhea | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 122 | EBV infection | Virus | 05 Viral Relief, 02 Core Support |
| 123 | Echinococcosis (hydatid) | Parasite | 14 Parasite Care, 02 Core Support |
| 124 | Echovirus | Virus | 05 Viral Relief, 02 Core Support |
| 125 | Endometrial cancer | Cancer | 15 Cancer Care, 02 Core Support |
| 126 | Endometriosis | Condition | 11 Reproductive Health, 02 Core Support |
| 127 | Enterobacter cloacae infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 128 | Enterobiasis (pinworm) | Parasite | 14 Parasite Care, 02 Core Support |
| 129 | Enterococcus faecalis infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 130 | Enterovirus (other) | Virus | 05 Viral Relief, 02 Core Support |
| 131 | Enterovirus 71 | Virus | 05 Viral Relief, 02 Core Support |
| 132 | Enterovirus D68 | Virus | 05 Viral Relief, 02 Core Support |
| 133 | Enterovirus infection | Virus | 05 Viral Relief, 02 Core Support |
| 134 | Epilepsy | Condition | 21 Brain Wave, 02 Core Support |
| 135 | Ewing sarcoma | Cancer | 15 Cancer Care, 02 Core Support |
| 136 | Eye strain / Computer vision syndrome | Condition | 49 Eyes Support, 10 Organ Health, 02 Core Support |
| 137 | Fascioliasis (rare) | Parasite | 14 Parasite Care, 02 Core Support |
| 138 | Fatty liver disease | Condition | 50 Liver Support, 02 Core Support |
| 139 | Fibromyalgia | Condition | 13 Arthritis/Pain, 22 Healing Mode |
| 140 | Filariasis hydrocele | Parasite | 14 Parasite Care, 02 Core Support |
| 141 | Fistula (general) | Condition | 16 Ulcer Care, 02 Core Support |
| 142 | Flea infestation | Parasite | 14 Parasite Care, 02 Core Support |
| 143 | Food poisoning (bacterial) | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |

| S/No | Disease/Pathogen | Type | Recommended mode(s) |
|-------------|---|-------------|--|
| 144 | Francisella tularensis (tularemia risk) | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 145 | Frozen shoulder | Condition | 13 Arthritis/Pain, 22 Healing Mode |
| 146 | Fungal otitis externa | Fungus | 04 Common Toxins/General Infection, 02 Core Support |
| 147 | Fungal sinusitis | Fungus | 04 Common Toxins/General Infection, 07 Dental Care |
| 148 | Fungal skin infection (unspecified) | Fungus | 04 Common Toxins/General Infection, 07 Dental Care |
| 149 | Fungal vaginitis | Fungus | 38 Vaginitis, 02 Core Support |
| 150 | Gallbladder cancer | Cancer | 15 Cancer Care, 02 Core Support |
| 151 | Gallbladder pain / gallstones support | Condition | 46 GallBldr Pain, 02 Core Support |
| 152 | Gastritis (bacterial) | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 153 | GERD | Condition | 10 Organ Health, 02 Core Support |
| 154 | Giardia chronic | Parasite | 14 Parasite Care, 02 Core Support |
| 155 | Giardiasis | Parasite | 14 Parasite Care, 02 Core Support |
| 156 | Gingivitis | Condition | 02 Core Support |
| 157 | GIST (gastrointestinal stromal tumor) | Cancer | 15 Cancer Care, 02 Core Support |
| 158 | Glaucoma | Condition | 49 Eyes Support, 10 Organ Health, 06 ENT Care, 02 Core Support |
| 159 | Glioblastoma | Cancer | 15 Cancer Care, 02 Core Support |
| 160 | Glioma (other) | Cancer | 15 Cancer Care, 02 Core Support |
| 161 | Gnathostomiasis | Parasite | 14 Parasite Care, 02 Core Support |
| 162 | Gonorrhea | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 163 | Gout | Condition | 13 Arthritis/Pain, 02 Core Support |
| 164 | Graves disease | Condition | 44 Thyroid Care, 25 Hormonal Care, 02 Core Support |
| 165 | Group B strep infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 166 | H. pylori infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 167 | Haemophilus ducreyi (chancroid) | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 168 | Haemophilus influenzae infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 169 | Haemorrhoids | Condition | 17 Piles/Gut, 02 Core Support |
| 170 | Hand foot mouth disease | Virus | 05 Viral Relief, 02 Core Support |
| 171 | Hangover | Condition | 50 Liver Support |
| 172 | Hantavirus infection | Virus | 05 Viral Relief, 02 Core Support |
| 173 | Hashimoto thyroiditis | Condition | 44 Thyroid Care, 25 Hormonal Care, 02 Core Support |

| S/No | Disease/Pathogen | Type | Recommended mode(s) |
|-------------|---|-------------|---|
| 174 | Headache (tension/cluster/sinus headache) | Condition | 47 Migraine Relief, 06 ENT Care (if sinus), 04 Common Toxins (if fever/cold), 02 Core Support |
| 175 | Heart failure | Condition | 08 Cardiac Care, 02 Core Support |
| 176 | Helicobacter-related ulcer | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 177 | Helminth infection (unspecified) | Parasite | 14 Parasite Care, 02 Core Support |
| 178 | Hepatitis A | Virus | 50 Liver Support, 05 Viral Relief, 02 Core Support |
| 179 | Hepatitis B | Virus | 50 Liver Support, 05 Viral Relief, 02 Core Support |
| 180 | Hepatitis C | Virus | 50 Liver Support, 05 Viral Relief, 02 Core Support |
| 181 | Hepatitis D | Virus | 50 Liver Support, 05 Viral Relief, 02 Core Support |
| 182 | Hepatitis E | Virus | 50 Liver Support, 05 Viral Relief, 02 Core Support |
| 183 | Hepatitis G/GBV-C | Virus | 50 Liver Support, 05 Viral Relief, 02 Core Support |
| 184 | Herpes simplex | Virus | 05 Viral Relief, 02 Core Support |
| 185 | Herpes zoster | Virus | 05 Viral Relief, 02 Core Support |
| 186 | HHV-6 | Virus | 05 Viral Relief, 02 Core Support |
| 187 | HHV-7 | Virus | 05 Viral Relief, 02 Core Support |
| 188 | HHV-8 (KSHV) | Virus | 05 Viral Relief, 02 Core Support |
| 189 | Histoplasmosis | Fungus | 04 Common Toxins/General Infection, 02 Core Support |
| 190 | HIV | Virus | 05 Viral Relief, 02 Core Support |
| 191 | Hodgkin lymphoma | Cancer | 15 Cancer Care, 02 Core Support |
| 192 | Hookworm disease | Parasite | 14 Parasite Care, 02 Core Support |
| 193 | Hospital-acquired pneumonia | Condition | 10 Organ Health, 04 Common Toxins/General Infection |
| 194 | HPV infection | Virus | 36 Cervix HPV (optional 11 Reproductive) |
| 195 | HTLV-1 | Virus | 05 Viral Relief, 02 Core Support |
| 196 | Human astrovirus | Virus | 05 Viral Relief, 02 Core Support |
| 197 | Human bocavirus | Virus | 05 Viral Relief, 02 Core Support |
| 198 | Human metapneumovirus | Virus | 05 Viral Relief, 02 Core Support |
| 199 | Human par echovirus | Virus | 05 Viral Relief, 02 Core Support |
| 200 | Hydatid disease | Parasite | 14 Parasite Care, 02 Core Support |
| 201 | Hymenolepiasis | Parasite | 14 Parasite Care, 02 Core Support |
| 202 | Hypertension | Condition | 08 Cardiac Care, 02 Core Support |
| 203 | Hypertensive urgency | Condition | 08 Cardiac Care, 02 Core Support |
| 204 | Hyperthyroidism | Condition | 44 Thyroid Care, 25 Hormonal Care, 02 Core Support |

| S/No | Disease/Pathogen | Type | Recommended mode(s) |
|-------------|---|-------------|--|
| 205 | Hypothyroidism | Condition | 44 Thyroid Care, 25 Hormonal Care, 02 Core Support |
| 206 | IBS | Condition | 17 Piles/Gut, 02 Core Support |
| 207 | Influenza A | Virus | 05 Viral Relief, 02 Core Support |
| 208 | Influenza B | Virus | 05 Viral Relief, 02 Core Support |
| 209 | Insomnia | Condition | 28 Sleep Relax |
| 210 | Interstitial lung disease | Condition | 10 Organ Health, 04 Common Toxins/General Infection |
| 211 | intestinal parasites (mixed) | Parasite | 14 Parasite Care, 02 Core Support |
| 212 | Intestinal worm infection (unspecified) | Parasite | 14 Parasite Care, 02 Core Support |
| 213 | Japanese encephalitis | Virus | 05 Viral Relief, 02 Core Support |
| 214 | Jaundice (non-obstructive) | Condition | 50 Liver Support, 02 Core Support |
| 215 | JC virus | Virus | 05 Viral Relief, 02 Core Support |
| 216 | Kala-azar (visceral leishmaniasis) | Parasite | 14 Parasite Care, 02 Core Support |
| 217 | Kidney stones | Condition | 09 Metabolic Care, 02 Core Support |
| 218 | Klebsiella granulomatis (donovanosis) | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 219 | Klebsiella pneumoniae infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 220 | Kyasanur Forest Disease | Virus | 05 Viral Relief, 02 Core Support |
| 221 | Laryngeal cancer | Cancer | 15 Cancer Care, 02 Core Support |
| 222 | Lassa fever (imported risk) | Virus | 05 Viral Relief, 02 Core Support |
| 223 | Legionella infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 224 | Legionellosis | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 225 | Leishmaniasis PKDL | Parasite | 14 Parasite Care, 02 Core Support |
| 226 | Leprosy | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 227 | Leptospira interrogans infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 228 | Leptospirosis | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 229 | Lipoma | Condition | 50 Liver Support, 02 Core Support |
| 230 | Listeria monocytogenes infection | Bacteria | 50 Liver Support, 04 Common Toxins, 02 Core Support |
| 231 | Liver Cirrhosis | Condition | 50 Liver Support, 02 Core Support |
| 232 | Liver Enlargement (Hepatomegaly) | Condition | 50 Liver Support, 09 Metabolic Care, 02 Core Support |
| 233 | Low back pain | Condition | 13 Arthritis/Pain, 02 Core Support |
| 234 | Lyme disease (rare/imported) | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 235 | Lymphatic filariasis | Parasite | 14 Parasite Care, 02 Core Support |
| 236 | Macular Degeneration (AMD) | Condition | 49 Eyes Support, 10 Organ Health, 02 Core Support |

| S/No | Disease/Pathogen | Type | Recommended mode(s) |
|-------------|------------------------------------|-------------|--|
| 237 | Malaria (P. falciparum) | Parasite | 14 Parasite Care, 02 Core Support |
| 238 | Malaria (P. vivax) | Parasite | 14 Parasite Care, 02 Core Support |
| 239 | Malaria relapse/fever (vivax) | Parasite | 14 Parasite Care, 02 Core Support |
| 240 | Measles | Virus | 05 Viral Relief, 02 Core Support |
| 241 | Melanoma | Cancer | 15 Cancer Care, 02 Core Support |
| 242 | Melioidosis | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 243 | Meningioma | Cancer | 15 Cancer Care, 02 Core Support |
| 244 | Menopause symptoms | Condition | 43 Menopause Help, 25 Hormonal Care |
| 245 | MERS-CoV (imported risk) | Virus | 05 Viral Relief, 02 Core Support |
| 246 | Metabolic syndrome | Condition | 02 Core Support |
| 247 | Microsporidiosis | Parasite | 14 Parasite Care, 02 Core Support |
| 248 | Migraine, Headache | Condition | 47 Migraine Relief, 22 Healing Mode, 02 Core Support |
| 249 | Morganella morganii infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 250 | Mpox | Virus | 05 Viral Relief, 02 Core Support |
| 251 | MRSA infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 252 | MSSA infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 253 | Mucormycosis | Fungus | 04 Common Toxins/General Infection, 07 Dental Care |
| 254 | Multidrug-resistant TB | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 255 | Multiple myeloma | Cancer | 15 Cancer Care, 02 Core Support |
| 256 | Mumps | Virus | 05 Viral Relief, 02 Core Support |
| 257 | Mycoplasma genitalium infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 258 | Mycoplasma pneumoniae infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 259 | Myiasis | Parasite | 14 Parasite Care, 02 Core Support |
| 260 | Naegleria fowleri infection | Parasite | 14 Parasite Care, 02 Core Support |
| 261 | Nasopharyngeal carcinoma | Cancer | 15 Cancer Care, 02 Core Support |
| 262 | Neisseria meningitidis infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 263 | Neuroblastoma | Cancer | 15 Cancer Care, 02 Core Support |
| 264 | Neurocysticercosis | Parasite | 14 Parasite Care, 02 Core Support |
| 265 | Nipah | Virus | 05 Viral Relief, 02 Core Support |
| 266 | Nocardia infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 267 | Non-Hodgkin lymphoma | Cancer | 15 Cancer Care, 02 Core Support |
| 268 | Non-typhoidal Salmonella infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |

| S/No | Disease/Pathogen | Type | Recommended mode(s) |
|-------------|----------------------------------|-------------|---|
| 269 | Norovirus gastroenteritis | Virus | 05 Viral Relief, 02 Core Support |
| 270 | Obesity | Condition | 02 Core Support |
| 271 | Onychomycosis | Fungus | 04 Common Toxins/General Infection, 07 Dental Care |
| 272 | Opisthorchiasis (rare) | Parasite | 14 Parasite Care, 02 Core Support |
| 273 | Oral candidiasis | Condition | 02 Core Support |
| 274 | Oral thrush | Fungus | 04 Common Toxins/General Infection, 07 Dental Care |
| 275 | Orientia tsutsugamushi infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 276 | Oropharyngeal cancer | Cancer | 15 Cancer Care, 02 Core Support |
| 277 | Osteoporosis | Condition | 45 Bone Health |
| 278 | Osteosarcoma | Cancer | 15 Cancer Care, 02 Core Support |
| 279 | Otitis externa | Condition | 02 Core Support |
| 280 | Otitis media (bacterial) | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 281 | Otitis media (recurrent) | Condition | 02 Core Support |
| 282 | Ovarian cysts | Condition | 42 Ovarian Cysts |
| 283 | Pancreatic cancer | Cancer | 15 Cancer Care, 02 Core Support |
| 284 | Parainfluenza | Virus | 05 Viral Relief, 02 Core Support |
| 285 | Parainfluenza type 1 | Virus | 05 Viral Relief, 02 Core Support |
| 286 | Parainfluenza type 2 | Virus | 05 Viral Relief, 02 Core Support |
| 287 | Parainfluenza type 3 | Virus | 05 Viral Relief, 02 Core Support |
| 288 | Parainfluenza type 4 | Virus | 05 Viral Relief, 02 Core Support |
| 289 | Parasites Flukes Liver | Parasite | 50 Liver Support, 14 Parasite Care, 02 Core Support |
| 290 | Paratyphoid fever | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 291 | Parvovirus B19 | Virus | 05 Viral Relief, 02 Core Support |
| 292 | PCOS | Condition | 11 Reproductive Health, 02 Core Support |
| 293 | Pediculosis (lice) | Parasite | 14 Parasite Care, 02 Core Support |
| 294 | Pemphigus vulgaris | Condition | 32 Pemphigus Vul, 12 Autoimmune, 02 Core Support |
| 295 | Penile cancer | Cancer | 15 Cancer Care, 02 Core Support |
| 296 | Peptic ulcer disease | Condition | 16 Ulcer Care, 02 Core Support |
| 297 | Periodontal disease | Condition | 02 Core Support |
| 298 | Periodontitis | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 299 | Peripheral artery disease | Condition | 08 Cardiac Care, 02 Core Support |
| 300 | Pertussis | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |

| S/No | Disease/Pathogen | Type | Recommended mode(s) |
|-------------|-----------------------------------|-------------|---|
| 301 | Pesticide Detox | Condition | 50 Liver Support, 23 Detox Body, 02 Core Support |
| 302 | Pharyngitis | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 303 | Pharyngitis (recurrent) | Condition | 02 Core Support |
| 304 | PID (Pelvic Inflammatory Disease) | Bacteria | 39 PID Adnexitis (optional 41 STIs RTIs) |
| 305 | Pimples | Bacteria | 34 Skin Wound, 04 Common Toxins |
| 306 | Pituitary tumor | Cancer | 15 Cancer Care, 02 Core Support |
| 307 | Plague | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 308 | Pneumococcal pneumonia | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 309 | Pneumocystis jirovecii pneumonia | Fungus | 04 Common Toxins/General Infection, 02 Core Support |
| 310 | Polio (wild/VDPV) | Virus | 05 Viral Relief, 02 Core Support |
| 311 | Prostatitis (chronic) | Condition | 11 Reproductive Health, 02 Core Support |
| 312 | Proteus mirabilis infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 313 | Protozoal diarrhea (unspecified) | Parasite | 14 Parasite Care, 02 Core Support |
| 314 | Protozoal infection (unspecified) | Parasite | 14 Parasite Care, 02 Core Support |
| 315 | Providencia infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 316 | Pseudomonas aeruginosa infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 317 | Pulmonary fibrosis | Condition | 13 Arthritis/Pain, 22 Healing Mode |
| 318 | Pyelonephritis | Bacteria | 33 UTI Kidney |
| 319 | Rabies | Virus | 05 Viral Relief, 02 Core Support |
| 320 | Recurrent UTI | Condition | 33 UTI Kidney |
| 321 | Renal cell carcinoma | Cancer | 15 Cancer Care, 02 Core Support |
| 322 | Retinoblastoma | Cancer | 15 Cancer Care, 02 Core Support |
| 323 | Rheumatic fever | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 324 | Rheumatoid arthritis | Condition | 13 Arthritis/Pain, 02 Core Support |
| 325 | Rhinovirus (other types) | Virus | 05 Viral Relief, 02 Core Support |
| 326 | Rhinovirus common cold | Virus | 05 Viral Relief, 02 Core Support |
| 327 | Rickettsia conorii infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 328 | Rickettsial fever (other) | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 329 | Rotavirus B | Virus | 05 Viral Relief, 02 Core Support |
| 330 | Rotavirus C | Virus | 05 Viral Relief, 02 Core Support |
| 331 | Rotavirus gastroenteritis | Virus | 05 Viral Relief, 02 Core Support |
| 332 | RSV | Virus | 05 Viral Relief, 02 Core Support |

| S/No | Disease/Pathogen | Type | Recommended mode(s) |
|-------------|--|-------------|---|
| 333 | RSV A | Virus | 05 Viral Relief, 02 Core Support |
| 334 | RSV B | Virus | 05 Viral Relief, 02 Core Support |
| 335 | Rubella | Virus | 05 Viral Relief, 02 Core Support |
| 336 | Sapovirus | Virus | 05 Viral Relief, 02 Core Support |
| 337 | SARS (risk) | Virus | 05 Viral Relief, 02 Core Support |
| 338 | Scabies | Parasite | 14 Parasite Care, 02 Core Support |
| 339 | Scarlet fever | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 340 | Schistosomiasis (rare) | Parasite | 14 Parasite Care, 02 Core Support |
| 341 | Sciatica | Condition | 13 Arthritis/Pain, 02 Core Support |
| 342 | Scrub typhus | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 343 | Seasonal coronavirus (unspecified) | Virus | 05 Viral Relief, 02 Core Support |
| 344 | Sepsis | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 345 | Sepsis syndrome | Condition | 02 Core Support |
| 346 | Serratia marcescens infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 347 | Shigella flexneri infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 348 | Shigella sonnei infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 349 | Shigellosis | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 350 | Sinusitis (bacterial) | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 351 | Soft tissue sarcoma | Cancer | 15 Cancer Care, 02 Core Support |
| 352 | Soil-transmitted helminth infection | Parasite | 14 Parasite Care, 02 Core Support |
| 353 | Squamous cell skin cancer | Cancer | 15 Cancer Care, 02 Core Support |
| 354 | Staph skin infection | Bacteria | 34 Skin Wound |
| 355 | Stenotrophomonas maltophilia infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 356 | Strep throat | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 357 | Stroke (hemorrhagic) | Condition | 18 Paralysis/Nerve, 02 Core Support |
| 358 | Stroke (ischemic) | Condition | 18 Paralysis/Nerve, 02 Core Support |
| 359 | Strongyloidiasis | Parasite | 14 Parasite Care, 02 Core Support |
| 360 | Syphilis | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 361 | Taenia solium cysticercosis | Parasite | 14 Parasite Care, 02 Core Support |
| 362 | Taeniasis | Parasite | 14 Parasite Care, 02 Core Support |
| 363 | Testicular cancer | Cancer | 15 Cancer Care, 02 Core Support |
| 364 | Tetanus | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |

| S/No | Disease/Pathogen | Type | Recommended mode(s) |
|-------------|--|-------------|---|
| 365 | Thymoma | Cancer | 15 Cancer Care, 02 Core Support |
| 366 | Thyroid cancer | Cancer | 15 Cancer Care, 02 Core Support |
| 367 | Tick infestation | Parasite | 14 Parasite Care, 02 Core Support |
| 368 | Tick-borne encephalitis (travel risk) | Virus | 05 Viral Relief, 02 Core Support |
| 369 | Tinea capitis | Fungus | 04 Common Toxins/General Infection, 02 Core Support |
| 370 | Tinea corporis | Fungus | 04 Common Toxins/General Infection, 02 Core Support |
| 371 | Tinea cruris | Fungus | 04 Common Toxins/General Infection, 02 Core Support |
| 372 | Tinea pedis (athlete foot) | Fungus | 04 Common Toxins/General Infection, 07 Dental Care |
| 373 | Tinea versicolor | Fungus | 04 Common Toxins/General Infection, 02 Core Support |
| 374 | Tonsillitis | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 375 | Toxocariasis | Parasite | 14 Parasite Care, 02 Core Support |
| 376 | Toxoplasmosis | Parasite | 14 Parasite Care, 02 Core Support |
| 377 | Toxoplasmosis (congenital) | Parasite | 14 Parasite Care, 02 Core Support |
| 378 | Tracheoesophageal fistula | Condition | 16 Ulcer Care, 02 Core Support |
| 379 | Trachoma | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 380 | Traveler diarrhea | Condition | 17 Piles/Gut, 02 Core Support |
| 381 | Treponema pallidum (syphilis) | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 382 | Trichinellosis (rare) | Parasite | 14 Parasite Care, 02 Core Support |
| 383 | Trichomoniasis | Parasite | 14 Parasite Care, 02 Core Support |
| 384 | Trichuriasis | Parasite | 14 Parasite Care, 02 Core Support |
| 385 | Tuberculosis | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 386 | Type 1 diabetes | Condition | 09 Metabolic Care, 02 Core Support |
| 387 | Type 2 diabetes | Condition | 09 Metabolic Care, 02 Core Support |
| 388 | Typhoid fever | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 389 | Ureaplasma urealyticum infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 390 | Urinary incontinence / enuresis / bedwetting | Condition | 40 Urinary Incont, 02 Core Support |
| 391 | Uterine fibroids | Condition | 11 Reproductive Health, 02 Core Support |
| 392 | UTI (E. coli) | Bacteria | 33 UTI Kidney (optional 04 Common Toxins) |
| 393 | Vaginal cancer | Cancer | 15 Cancer Care, 02 Core Support |
| 394 | Vaginitis / vaginosis | Condition | 38 Vaginitis, 02 Core Support |
| 395 | Varicella (chickenpox) | Virus | 05 Viral Relief, 02 Core Support |

| S/No | Disease/Pathogen | Type | Recommended mode(s) |
|-------------|-------------------------------------|-------------|---|
| 396 | Varicose veins | Condition | 22 Healing Mode, 02 Core Support |
| 397 | Venous insufficiency | Condition | 08 Cardiac Care, 02 Core Support |
| 398 | Ventilator-associated pneumonia | Condition | 10 Organ Health, 04 Common Toxins/General Infection |
| 399 | Vibrio parahaemolyticus infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 400 | Vibrio vulnificus infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 401 | Viral bronchiolitis | Virus | 05 Viral Relief, 02 Core Support |
| 402 | Viral conjunctivitis | Virus | 05 Viral Relief, 02 Core Support, 49 Eyes Support |
| 403 | Viral encephalitis (unspecified) | Virus | 05 Viral Relief, 02 Core Support |
| 404 | Viral exanthem (unspecified) | Virus | 05 Viral Relief, 02 Core Support |
| 405 | Viral fever (unspecified) | Virus | 05 Viral Relief, 02 Core Support |
| 406 | Viral gastroenteritis (unspecified) | Virus | 05 Viral Relief, 02 Core Support |
| 407 | Viral meningitis (enteroviral) | Virus | 05 Viral Relief, 02 Core Support |
| 408 | Viral myocarditis (enteroviral) | Virus | 05 Viral Relief, 02 Core Support |
| 409 | Vitamin D deficiency | Condition | 02 Core Support |
| 410 | Vitiligo | Condition | 50 Liver Support, 12 Autoimmune, 02 Core Support |
| 411 | VRE (Enterococcus faecium) | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 412 | Vulvar cancer | Cancer | 15 Cancer Care, 02 Core Support |
| 413 | West Nile fever | Virus | 05 Viral Relief, 02 Core Support |
| 414 | Wilms tumor | Cancer | 15 Cancer Care, 02 Core Support |
| 415 | Yellow Fever | Virus | 50 Liver Support, 05 Viral Relief, 02 Core Support |
| 416 | Yersinia enterocolitica infection | Bacteria | 04 Common Toxins/General Infection, 02 Core Support |
| 417 | Zika | Virus | 05 Viral Relief, 02 Core Support |

Product Warranty (Hardware)

Limited Hardware Warranty

The manufacturer warrants that the Zenthera Zapper hardware (circuit board, LCD, buttons, enclosure, and buzzer) will be free from defects in materials and workmanship under normal use for a period of **one (1) year** from the date of purchase.

The power adapter, cables, and electrodes are covered by a limited warranty for **three (3) months** from the date of purchase.

Conductive gel is a consumable item and is provided **without warranty**. Once consumed, it can be purchased separately from nearby medical supply shops.

This warranty covers:

- Manufacturing defects in electronic components
- LCD display failures not caused by physical damage
- Button malfunction due to manufacturing defects
- Power circuit failures under normal operating conditions

This warranty does NOT cover:

- Damage from misuse, abuse, or improper handling
- Water damage or exposure to moisture beyond specified IP rating
- Battery deterioration (batteries are consumable items)
- Normal wear and tear
- Damage from unauthorized modifications or repairs
- Use with incorrect voltage or current specifications

Warranty Claim Process:

Contact the volunteer or our website www.thetempleofhealing.org or write a email on CureYourSelfIndia@gmail.com with description of defect. We will coordinate with you and will repair or replace defective units at its sole discretion.

Disclaimer of Medical Warranty

NO MEDICAL WARRANTY

This device is provided "AS IS" for experimental and research purposes only. The manufacturer makes **NO WARRANTY, EXPRESS OR IMPLIED**, regarding:

- The device's effectiveness for any health condition
- The safety or efficacy of frequency-based therapy
- Medical outcomes or health improvements
- Accuracy of frequency programs for any stated purpose

The manufacturer explicitly disclaims:

- Any warranty of merchantability for medical use
- Any warranty of fitness for a particular medical purpose
- Any claims that this device can diagnose, treat, cure, or prevent disease

This device is **NOT** a medical device. It has **NOT** been evaluated or approved by any regulatory health authority (DCGI, CDSCO, NRA, FDA, CE, TGA, or equivalent).

Limitation of Liability

LIMITATION OF LIABILITY

To the maximum extent permitted by law, the manufacturer's liability is limited to:

- The purchase price of the device
- Repair or replacement of defective hardware (within warranty period)

The manufacturer shall **NOT** be liable for:

- Personal injury or health complications arising from device use
- Consequential, incidental, special, or punitive damages
- Medical expenses or healthcare costs
- Loss of time, inconvenience, or loss of use
- Reliance on frequency programs for health decisions
- Interference with medical devices (including pacemakers)

USER ASSUMES ALL RISK associated with using this experimental device.

Safety Disclaimer (Critical)

SAFETY WARNING & DISCLAIMER

 **DO NOT USE if you have:**

- A pacemaker or any implanted electronic medical device
- Are pregnant or suspect pregnancy

- Have epilepsy or seizure disorders
- Have a heart condition without physician clearance

This device generates electromagnetic frequencies. The long-term health effects are unknown. Use is entirely at your own risk.

The manufacturer provides this device for educational and experimental purposes. Users are solely responsible for their health decisions and should consult qualified healthcare professionals before using frequency-based therapies.

Software/Firmware Disclaimer

FIRMWARE & FREQUENCY PROGRAMS DISCLAIMER

The frequency programs stored in this device are compiled from various public-domain sources and historical Rife frequency lists. The manufacturer makes **NO REPRESENTATION** regarding:

- The accuracy of frequency values
- The safety of any specific frequency
- The therapeutic value of any program
- The scientific validity of frequency-based approaches

Frequency programs may be updated via firmware. Users accept all risks associated with firmware updates and frequency program changes.

No Endorsement Disclaimer

NO MEDICAL ENDORSEMENT

References to health conditions, symptoms, or organs in mode labels (e.g., "Cancer Support," "Cardiac Care," "UTI Kidney") are:

- Descriptive labels for frequency list organization only
- NOT medical claims or treatment promises
- NOT endorsed by medical professionals
- NOT validated by clinical studies

These labels reflect historical Rife frequency list categories and do NOT constitute medical advice or treatment recommendations.

<!-- End of File -->